

Dental Equipment Retail Price List

April 2026_{v2}

See our new website
belmontdental.co.uk

Dental Equipment Retail Price List April 2026 CONTENTS

- 2 **FREE 5 Year Warranty**
- 4 **New Optional Implantology Motor System**
- 6 **Eurus S1**
Treatment Centres / Flexible Options / Factory Built-in Options
- 12 **Eurus S4**
Treatment Centres / Factory Built-in Options
- 16 **Eurus S6**
Treatment Centres / Flexible Options / Factory Built-in Options
- 22 **Eurus S8**
Treatment Centres / Factory Built-in Options
- 30 **Voyager III Plus**
Treatment Centres / Flexible Options / Factory Built-in Options
- 34 **Flexible Delivery Options**
Voyager 'M' and 'MC' / Eurus 'M' and 'MC'
- 38 **Dental Chairs**
Eurus Swivel / Eurus Contour
- 40 **Intra Oral X-rays**
097 Belray II / Belmont Touch Phot-X IIS DC LCD
- 42 **Dental Operating Lights**
100 LED Series / 900 LED Series
- 44 **Dental Stools**
For Operators & Assistants
- 46 **Upholstery**
Standard / Ultrasoft / Pearlescent
- 50 **New Dental Educational Treatment Unit**
- 51 **Hygiene & Care**
Cleaning and maintenance range
- 52 **Service Kits**
For Treatment Centres, Modules & Mobile Carts

Belmont

All prices listed are
shown exclusive of VAT.

The manufacturers reserve the right to alter
the prices and specifications without prior notice.

See our
new website:



Tell me about the
Eurus touch-screen:



See brochure for
Dental Educational Unit:



FREE

5 Year Warranty

All current models of Belmont dental chairs have been through **vigorous durability testing** prior to mass production ensuring total reliability for the operator.

Each chair design has been tested to a ten-year life cycle of use, with 108,000 movements fully loaded to the chair maximum safe working patient weight to regulation ISO 7494-1 specification.

(Which is equivalent to: 10 Years x 240 Days per year x 15 Patients per day x 3 chair movements per patient).



Takara Belmont (UK) Ltd are registered with RECYCLINGLIVES SERVICES for:
Electrical WEEE Waste
Registration No: WEE/DH0151QX
and Battery and Accumulator Waste
Registration No: BPRN12034

Takara Belmont (UK) Ltd Warranty Policy

5 Year Warranted Products

All Belmont treatment centres, individual chairs, units and dental lights are warranted for a period of five years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions.

Upholstery, paintwork, plating and all handpiece hoses are not warranted, unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery. Failure due to natural wear and tear of any component is not warranted.

OEM factory build-in options (ie: fibre optics, scalers, motors, separators, CVS Kits etc) carry a warranty as recommended by the individual manufacturer.

2 Year Warranted Products

All Belmont intra oral x-rays are warranted for two years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions. Paintwork and plating are not warranted unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery.

Takara Belmont (UK) Ltd are not responsible for any failure caused by misuse of the equipment or damage caused due to poor installation by a non-authorised company.

The warranty policy applies to equipment purchased from and installed by an authorised Takara Belmont (UK) Ltd dealer. Labour will be charged on warranty work after the first year.

Hygiene & Care

In order to maintain & prolong the longevity of your Belmont equipment, we highly recommend using our Hygiene & Care products. (See page 51)

PLEASE NOTE: ISO 13485 Certification

All of our dental products are supplied with a full EU Declaration of Conformity and certification to the standard complying with ISO 13485. Any modifications or alterations made after certification may affect your warranty as the responsibility for EU Declaration of Conformity and electrical safety will transfer to the modifier.

Belmont Dental are delighted to bring implantology integration as a built-in option to **all Eurus treatment centres**

IN PARTNERSHIP WITH BIEN-AIR

Implantology Built-In System

Optional for the whole range of E-type Eurus Units and all operational via the intuitive touch-screen

Bien-Air MX-i LED 3rd Generation Micro Motor

CA 20:1 Contra Angle Micro Series Handpiece

Peristaltic Irrigation Pump

Motor Stand

Irrigation Line Stand

Irrigation Lines (Pack of 10)

Irrigation Clip

Foot Stalk for Program Control
(assembles onto the unit foot control)



Important Note: This option must be purchased with our Eurus products from the factory. This cannot be retro fitted on site.

TREATMENT CENTRES (Over-the-Patient delivery system)	With Unit Mounted D-921 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S1 'Surgery System' Holder Type (A)	£22,805	£23,280	£21,005
Eurus S1 'Surgery System' Holder Type (E)	£25,915	£26,390	£24,115
Eurus S1 'Surgery System' CM Rod Type (A)	£24,090	£24,565	£22,290
Eurus S1 'Surgery System' CM Rod Type (E)	£27,210	£27,685	£25,410
Eurus S1 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£23,375	£23,850	£21,575
Eurus S1 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£26,485	£26,960	£24,685
When light sold with package above the saving over recommended retail price is:	Save £1,295	Save £1,020	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Chair Mounted Over-Patient Holder & Rod Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation

consists of ~

Chair section features include:

The Eurus S1 has a maximum safe working patient weight of 200kg (ISO 7494-1)
Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Integrated service centre
Synchronised axis backrest & seat movement
Integrated foot stalks for chair control
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Autoclavable unit handle cover, handpiece holders
unit tray silicon mat, handpiece rest (Rod Type)
Easy-clean membrane switch panel
Doctor handpiece control touch screen with water low warning indicator
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Electronic foot control for accurate turbine speed control and integrated lifting handle

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with 2 fibre-optic outlets
Optional upgrade to WaveOne endodontic system
Optional upgrade to Implantology Motor System

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation
921 LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Bottled water system with water low warning
LED on cup filler spout casing
Durr spittoon valve



Eurus S1 over-the-patient holder type

Eurus S1 Treatment Centres



Flexible Treatment Centres

Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation

consists of ~

Chair section features include:

The Eurus S1 has a maximum safe working patient weight of 200kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Swing-out armrest to 90 or 135 degrees

Integrated service centre

Synchronised axis backrest & seat movement

Integrated foot stalks for chair control

Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurus Cabinet Module or Mobile Cart delivery systems A & E (please specify)

'A' version features include:

Cabinet mounting bracket and arm or height adjustable mobile cart base

Autoclavable unit handle cover, handpiece holders and unit tray silicon mat

Easy clean membrane switch panel

Doctor handpiece control touch screen

3-way autoclavable syringe

3 Fibre optic & 1 non optic handpiece outlets with individual air & water controls

First priority handpiece selection

Electronic foot control for accurate turbine speed control and integrated lifting handle

Handpiece flush-out

Handpiece oil mist collector

Bottled water system

'E' version additional features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with 2 fibre-optic outlets

Optional upgrade to WaveOne endodontic system

Optional upgrade to Implantology Motor System

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation

921 LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Bottled water system with water low warning

LED on cup filler spout casing

Durr spittoon valve

Important:

Maximum loom length between module or cart to the chair is 6.7metres (standard is 2.7metres).

Please read before placing your order:

When ordering a Flexible Treatment Centre, it is imperative that you order the correct item – either a **Cabinet Mounted Module or Mobile Cart**. If the **incorrect item** has been entered onto your purchase order and delivered to you, the cuspidor and delivery unit cannot be modified and **Takara Belmont (UK) Ltd will not accept returns**. Two completely new items will need to be manufactured to the correct specification, which will cause delays for your installation at the dental practice, and cost implications may be incurred.

Eurus S1



Eurus S1 shown with seamless Upholstery

Eurus S1 Factory Built-in Options



for Eurus S1	
£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,505	Additional Bien Air MX2 fibre-optic micro motor *
£3,930	Upgrade to Bien Air Implantology Motor System
£1,130	Upgrade to Bien Air WaveOne System with Contra Angle
£1,680	Additional NSK NLX-Plus fibre-optic micro motor *
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic scaler
£1,830	NSK Varios 170 optic scaler
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£3,070	Durr CAS 1 Combi-Separator with rinse
£1,730	Acteon C20 HD Intra-Oral Camera
£995	Wireless unit foot controller (when factory fitted)
£750	Monitor mount
£435	Additional 3-way autoclavable syringe (assistant side)
£210	Additional easily detachable HV ejector hose (including ball joint for greater ease of operation)
£195	Large tray (double size for Rod Type)
£100	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£250	Flushing Tray A (Holder Type)
£250	Flushing Tray B (Rod Type)
£280	Annual service Kit 7

* Additional micro motors must be the same model as the primary micro motor.



New Flushing Tray – suitable for all treatment centres

Eurus S4

TREATMENT CENTRES (Over-the-Patient delivery system)	With Unit Mounted D-921 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S4 'Surgery System' Pedestal Holder Type (A)	£26,170	£26,645	£24,370
Eurus S4 'Surgery System' Pedestal Holder Type (E)	£28,600	£28,870	£26,595
Eurus S4 'Surgery System' Pedestal Rod Type (A)	£27,095	£27,570	£25,295
Eurus S4 'Surgery System' Pedestal Rod Type (E)	£30,560	£30,830	£28,555
When light sold with package above the saving over recommended retail price is:	Save £1,295	Save £1,020	£0

Please state light model number when ordering treatment centre

Available to order with power headrest

Please refer to p.40 for unit mounted Intra Oral X-ray

Please see swatch brochure for upholstery options

Eurus S4 Treatment Centres

with folding and extending legrest



Pedestal Mounted Over-Patient Holder & Rod Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation

consists of ~

Chair section features include:

The Eurus S4 has a folding and extending legrest with a maximum safe working patient weight of 200kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Integrated foot stalks for chair control

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Autoclavable unit handle cover, handpiece holders

unit tray silicon mat, handpiece rest (Rod Type)

Easy-clean membrane switch panel

Doctor handpiece control touch screen with water low warning indicator

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate turbine speed control and integrated lifting handle

Integrated service centre

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with

auto-reverse, Endo control and LCD display, with 2 fibre-optic outlets

Optional upgrade to WaveOne endodontic system

Optional upgrade to Implantology Motor System

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation

921 LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Bottled water system with water low warning

LED on cup filler spout casing

Durr spittoon valve

Eurus S4



Eurus S4 highlighting pedestal and folding legrest chair

Eurus S4 Factory Built-in Options



for Eurus S4

£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,505	Additional Bien Air MX2 fibre-optic micro motor *
£3,930	Upgrade to Bien Air Implantology Motor System
£1,130	Upgrade to Bien Air WaveOne System with Contra Angle
£1,680	Additional NSK NLX-Plus fibre-optic micro motor *
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic scaler
£1,830	NSK Varios 170 optic scaler
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£3,070	Durr CAS 1 Combi-Separator with rinse
£1,730	Acteon C20 HD Intra-Oral Camera
£750	Monitor mount
£435	Additional 3-way autoclavable syringe (assistant side)
£210	Additional easily detachable HV ejector hose (including ball joint for greater ease of operation)
£195	Large tray (double size for Rod Type)
£100	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£250	Flushing Tray A (Holder Type)
£250	Flushing Tray B (Rod Type)
£280	Annual service Kit 7

* Additional micro motors must be the same model as the primary micro motor.



New Flushing Tray –
suitable for all treatment centres

TREATMENT CENTRES (Below-the-Patient delivery system)	With Unit Mounted D-921 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S6 'Surgery System' Holder Type (A)	£28,830	£29,305	£27,030
Eurus S6 'Surgery System' Holder Type (E)	£31,985	£32,460	£30,185
Eurus S6 'Vac-Pac System' Holder Type (A)	£23,955	£24,430	£22,155
Eurus S6 'Vac-Pac System' Holder Type (E)	£27,075	£27,550	£25,275
Eurus S6 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£28,150	£28,625	£26,350
Eurus S6 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£31,230	£31,705	£29,430
When light sold with package above the saving over recommended retail price is:	Save £1,295	Save £1,020	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Eurus S6 Treatment Centres

with folding and extending legrest



Surgery System (A and E)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg (ISO 7494-1) (Vac-Pac 200kg)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Taller backrest cushion with additional lumbar support for greater patient comfort

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray with pneumatic release and lock button

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece

holders unit tray silicon mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets

with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate

instrument operation with integrated lifting handle

Bottle water system

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with

auto-reverse, Endo control and LCD display, with 2 fibre-optic outlets

Optional upgrade to WaveOne endodontic system

Optional upgrade to Implantology Motor System

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with

90 degree rotation

Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left side

921 LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the cup filler spout and touch-screen

Durr spittoon valve

'Vac-Pac System' (A and E)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

Assistant's console

Easy-access detachable filter

Autoclavable HVE & saliva ejector

Funnel connection hose



Eurus S6 Treatment Centres

with folding and extending legrest



Flexible Treatment Centres Type (A and E)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Chair armrest folds away 180 degrees towards the back of the chair for easy infirm patient entry
Synchronised axis backrest & seat movement
Taller backrest cushion with additional lumbar support for greater patient comfort
Integrated foot stalks for chair control
Integrated service centre
Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurus Cabinet Module or Mobile Cart delivery systems A & E (please specify)

'A' version features include:

Cabinet mounting bracket and arm or height adjustable mobile cart base
Autoclavable unit handle cover, handpiece holders and unit tray silicon mat
Easy clean membrane switch panel
Doctor handpiece control touch screen
3-way autoclavable syringe
3 Fibre optic & 1 non optic handpiece outlets with individual air & water controls
First priority handpiece selection
Wireless electronic foot control for accurate instrument operation with integrated lifting handle
Handpiece flush-out
Handpiece oil mist collector
Bottled water system

'E' version additional features include:

Bien-Air MX2 fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with 2 fibre-optic outlets
Optional upgrade to WaveOne endodontic system
Optional upgrade to Implantology Motor System

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 90 degree rotation
Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left side
921 LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
2 x HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Chair movement lock indicators located on the cup filler spout and touch-screen
Durr spittoon valve

Important:

Maximum loom length between module or cart to the chair is 6.7metres (standard is 2.7metres).

Please read before placing your order:

When ordering a Flexible Treatment Centre, it is imperative that you order the correct item – either a **Cabinet Mounted Module or Mobile Cart**. If the **incorrect item** has been entered onto your purchase order and delivered to you, the cuspidor and delivery unit cannot be modified and **Takara Belmont (UK) Ltd will not accept returns**. Two completely new items will need to be manufactured to the correct specification, which will cause delays for your installation at the dental practice, and cost implications may be incurred.

Eurus S6



Eurus S6 shown with Luxury Upholstery

Eurus S6 Factory Built-in Options



for Eurus S6

£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,505	Additional Bien Air MX2 fibre-optic micro motor *
£3,930	Upgrade to Bien Air Implantology Motor System
£1,130	Upgrade to Bien Air WaveOne System with Contra Angle
£1,680	Additional NSK NLX-Plus fibre-optic micro motor *
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic scaler
£1,830	NSK Varios 170 optic scaler
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£3,070	Durr CAS 1 Combi-Separator with rinse
£1,730	Acteon C20 HD Intra-Oral Camera
£750	Monitor mount
£435	Additional 3-way autoclavable syringe (assistant side)
£100	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£250	Flushing Tray A (Holder Type)
£280	Annual Service Kit 7

* Additional micro motors must be the same model as the primary micro motor.

Eurus S8

Base mounted ambidextrous Holder Type

TREATMENT CENTRES (Below-the-Patient base mounted delivery system)	With Unit Mounted D-921LR Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S8 'Surgery System' Holder Type (A)	£24,930	£25,405	£23,130
Eurus S8 'Surgery System' Holder Type (E)	£28,035	£28,510	£26,235
Eurus S8 'Vac-Pac System' Holder Type (A)	£22,550	£23,025	£20,750
Eurus S8 'Vac-Pac System' Holder Type (E)	£25,660	£26,135	£23,860
When light sold with package above the saving over recommended retail price is:	Save £1,295	Save £1,020	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Eurus S8 Treatment Centres

base mounted delivery system



Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) &

auto-return plus 6 operator configurations

Swing-out armrest to 90 or 135 degrees

Synchronised axis backrest & seat movement

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray with 3 settings (Holder Type)

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders & tray silicon mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets

with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate

instrument operation with integrated lifting handle

plus chair manual control movements (Holder Type)

Bottle water system

Unit section (E) Doctor side additional

features include:

NSK NLX-Plus fibre optic micro motor

Pre-set micro motor & torque control with

auto-reverse & touch-screen display

2 fibre-optic outlets

Implantology and WaveOne Systems are available as an additional option when the Bien Air MX2 micro motor is fitted

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with

90 degree rotation

Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left and right sides

921LR LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the

cup filler spout and touch-screen

Durr spittoon valve

Eurus S8

Ambidextrous Holder & Rod Type Swing

TREATMENT CENTRES (Chair mounted Swing Type delivery system)	With Unit Mounted D-921LR Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S8 'Surgery System' CM Holder Type Swing (A)	£24,930	£25,405	£23,130
Eurus S8 'Surgery System' CM Holder Type Swing (E)	£28,035	£28,510	£26,235
Eurus S8 'Vac-Pac System' CM Holder Type Swing (A)	£22,550	£23,025	£20,750
Eurus S8 'Vac-Pac System' CM Holder Type Swing (E)	£25,660	£26,135	£23,860
Eurus S8 'Surgery System' CM Rod Swing (A)	£27,570	£28,045	£25,770
Eurus S8 'Surgery System' CM Rod Swing (E)	£30,595	£31,070	£28,795
Eurus S8 'Vac-Pac System' CM Rod Swing (A)	£24,525	£25,000	£22,725
Eurus S8 'Vac-Pac System' CM Rod Swing (E)	£27,555	£28,030	£25,755
When light sold with package above the saving over recommended retail price is:	Save £1,295	Save £1,020	£0

For further details please see page 61 in our Dental Equipment brochure for Eurus CM Holder Type Swing

Swing Type Vac Pac is identical to base mounted unit on assistants side (see page 26)

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Eurus S8 Treatment Centres

chair mounted swing type delivery system



Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Swing-out armrest to 90 or 135 degrees

Synchronised axis backrest & seat movement

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders

& tray silicon mat / Rod Type handpiece place mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets

with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate

instrument operation with integrated lifting handle

Bottle water system

Unit section (E) Doctor side additional

features include:

NSK NLX-Plus fibre optic micro motor

Pre-set micro motor & torque control with

auto-reverse & touch-screen display

2 fibre-optic outlets

Implantology and WaveOne Systems are available as an additional option when the Bien Air MX2 micro motor is fitted

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with

90 degree rotation

Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left and right sides

921LR LED dental light with detachable

autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the

cup filler spout and touch-screen

Durr spittoon valve

The CM Holder Type Swing treatment centres are identical in specification to the standard S8, and have a chair-mounted swing around delivery arm instead of the base mounted arm type

Eurus S8

Eurus S8 Vac-Pac Treatment Centre



Eurus S8 Vac-Pac shown with Seamless Upholstery

Eurus S8 Vac-Pac Treatment Centres



Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg (ISO 7494-1)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Swing-out armrest to 90 or 135 degrees

Synchronised axis backrest & seat movement

Integrated foot stalks for chair control (Swing Type)

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders & tray silicon mat / Rod Type handpiece place mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate instrument operation with integrated lifting handle plus chair manual control movements (Base Mounted Type)

Bottle water system

Unit section (E) Doctor side additional features include:

NSK NLX-Plus fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse & touch-screen display

2 fibre-optic outlets

Implantology and WaveOne Systems are available as an additional option when the Bien Air MX2 micro motor is fitted

All Eurus S8 Vac-Pac Systems (A and E)

(Ambidextrous)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-access detachable filter

Autoclavable HVE & saliva ejector

Funnel connection hose

Chair control auto & manual touch pad

Assistant tray

Eurus S8

Eurus S8 Holder Type
treatment centre



Eurus S8 Factory Built-in Options



for Eurus S8	
£0	Bien Air MX2 micro motor (in place of NSK NLX-Plus micro motor on 'E' version only)
£1,505	Additional Bien Air MX2 fibre-optic micro motor *
£3,930	Upgrade to Bien Air Implantology Motor System
£1,130	Upgrade to Bien Air WaveOne System with Contra Angle
£1,680	Additional NSK NLX-Plus fibre-optic micro motor *
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic scaler
£1,830	NSK Varios 170 optic scaler
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£3,070	Durr CAS 1 Combi-Separator with rinse
£1,730	Acteon C20 HD Intra-Oral Camera
£360	Monitor mount
£435	Additional 3-way autoclavable syringe (assistant side)
£195	Large tray (double size for Rod Type)
£100	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£250	Flushing Tray A (Holder Type)
£250	Flushing Tray B (Rod Type)
£280	Annual Service Kit 7

* Additional micro motors must be the same model as the primary micro motor.

Voyager III Plus Ambidextrous

TREATMENT CENTRES Below-the-Patient delivery system and flexible options	With Unit Mounted D-119 Light	With Ceiling Mounted D-112 Light	Without Dental Light
Voyager III Plus 'Surgery System'	£16,130	£16,550	£14,720
Voyager III Plus 'SP' Chair with ambidextrous cuspidor with Voyager Cabinet Mounted Module or Mobile Cart delivery system	£16,130	£16,550	£14,720
When light sold with package above the saving over recommended retail price is:	Save £895	Save £910	£0
Voyager III Plus 'SP' Chair with ambidextrous cuspidor (without operating light)			£11,380

Please state light model number when ordering treatment centre

MODULE & MOBILE CART OPTIONS

Voyager 'M' & 'MC' Module and Mobile Cart delivery system	£3,340
--	--------

Voyager III Plus Treatment Centres & Flexible Systems



'Surgery System'

(Ambidextrous)

with D-119LR unit-mounted Operating Light

consists of ~

Chair section features include:

The Voyager III Plus has a maximum safe working patient weight of 135kg (ISO 7494-1)

New twin-articulating headrest system

2 preset positions, Last Position Memory (LP)

& Auto-Return

Swing-out armrests (left & right)

Integral service centre with air & water filter regulators

8-way foot control

Base-mounted swing-around delivery

system features include:

3 handpiece outlets with individual air & water control

3-way autoclavable syringe

Table safety lock device

Quick flush-out system

Bottled water system

Stainless steel instrument tray

Non-slip mat

Chair-mounted rotatable cuspidor unit

features include:

Porcelain cuspidor bowl

Manual cup-filler & bowl-flush

1 detachable high volume evacuator

1 detachable saliva ejector

Easily-accessed & detached solids collector

Durr spittoon valve

'SP'

(Ambidextrous)

Voyager III Plus Dental Chair with cuspidor

consists of ~

Chair section features include:

The Voyager III Plus has a maximum safe working patient weight of 135kg (ISO 7494-1)

New twin-articulating headrest system

2 preset positions, Last Position Memory (LP)

& Auto-Return

Swing-out armrests (left & right)

Integral service centre with air & water filter regulators

8-way foot control

Chair-mounted rotatable cuspidor unit

features include:

Porcelain cuspidor bowl

Manual cup-filler & bowl-flush

1 detachable high volume evacuator

1 detachable saliva ejector

Easily-accessed & detached solids collector

Durr spittoon valve

for Voyager

£2,635	Bien Air Optima MCX LED (inc. standard fibre-optic micro motor, turbine hose & touch pad control)
£2,635	NSK NLX-Nano fibre-optic pack (inc. NLX-Nano fibre-optic micro motor, turbine hose & touch pad control)
£3,210	NSK NLX-Plus multi-pad endo pack (inc. NLX-Plus fibre-optic micro motor, turbine hose & multi-pad Endo control)
£870	Bien Air fibre-optic kit (inc. 1 turbine hose)
£315	Additional fibre-optic outlet
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic basic scaler
£1,830	NSK Varios 170 optic scaler (with touch pad control)
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£435	Additional 3-way autoclavable syringe on Assistant's side inc. bottled water system
£0	Standard (10x10) Module Arm (alternative to Voyager Arm System)
£545	Touch pad assembly for chair control fitted to Operator's table
£765	Monitor mount
£1,730	Acteon C20 HD Intra-Oral Camera
£95	Extra-long umbilical (up to 4 metres for Mobile Cart only)
£100	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£130	Floor box fitted to cart base
£250	Flushing Tray A (Holder Type)
£270	Annual Service Kit 1 (for Voyager I, II, III & III Plus Treatment Centres)
£175	Annual Service Kit 2 (for Voyager Module and Cart)

* Additional micro motors must be the same model as the primary micro motor.

VOYAGER III PLUS

Find out how we've made a **30 year**
best-seller **even better**



Belmont

Perfecting the art of dentistry

belmontdental.co.uk

Flexible Delivery Options

MODULE & MOBILE CART OPTIONS

Voyager 'M' & 'MC' Module and Mobile Cart delivery system	£3,340
Eurus 'M' & 'MC' (A) Module and Mobile Cart delivery system	£6,510
Eurus 'M' & 'MC' (E) Module and Mobile Cart delivery system	£9,620

Voyager 'M' & 'MC' **Factory Built-in Options**
Please refer to page 32 for prices

Eurus 'M' & 'MC' **Factory Built-in Options**
Please refer to page 37 for prices

Universally compatible
Flexible Options



Voyager 'M' & 'MC'

Module delivery system features include:

Mounting arms & bracket

Mobile Cart delivery system with
height-adjustable mobile base
features include:

3 handpiece outlets with individual air & water control

3-way autoclavable syringe

Bottled water system

Stainless steel instrument tray

Non-slip mat

Service centre with air filter regulator

Eurus 'M' & 'MC'

Eurus Cabinet Mounted Module and Mobile Cart Standalone or with Eurus Chair options

Unit (A) Features include:

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders, unit tray silicon mat.

Easy-clean membrane switch panel for chair control
3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Electronic foot control (wired) for accurate turbine speed control and integrated lifting handle

Bottled water system

Additional chair connecting loom length of 4m will be supplied when ordered with a Eurus Chair

Unit (E) with additional options.

Additional features & new integrated options include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display and two fibre optics

Available optional up-grade treatment functions:

WaveOne endodontic system

Implantology complete system including, Bien-Air MX-i micro motor and CA 20:1 handpiece, irrigation pump and stand

OR

Unit section (E)

Additional features include:

NSK NLX-Plus fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse & touch-screen display

2 fibre-optic outlets

Eurus Mobile Cart and Swivel Chair



for Eurus Module & Cart

£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,505	Additional Bien Air MX2 fibre-optic micro motor *
£3,930	Upgrade to Bien Air Implantology Motor System
£1,130	Upgrade to Bien Air WaveOne System with Contra Angle
£1,680	Additional NSK NLX-Plus fibre-optic micro motor *
£2,260	Cavitron G-139 scaler
£1,545	NSK Varios 170 non-optic scaler
£1,830	NSK Varios 170 optic scaler
£1,625	Satelec Newtron scaler
£1,915	Satelec Newtron LED scaler
£1,730	Acteon C20 HD Intra-Oral Camera
£995	Wireless unit foot controller (when factory fitted)
£25	Connecting lead (Comm.) for Eurus Chair to Cart/Cabinet Mount (5m) - 1A0WLAA0
£280	Annual service Kit 7

** Additional micro motors must be the same model as the primary micro motor.*

*Eurus Cabinet Mounted Module and Mobile Cart
Standalone or with Eurus Chair options*



Dental Chairs

WITH FOLDING AND EXTENDING LEGREST

Eurus Swivel Chair

Standalone chair (or linked to Eurus Cabinet Mounted Module or Cart) *with* Twin-articulating headrest & cable free base rotation 105 degrees (L & R), available in all seamless and luxury upholstery

£9,770

Eurus Swivel PH Chair

with Power headrest (special order)
available in all seamless and luxury upholstery

£11,125

WITH FIXED LEGREST

Eurus Contour Chair

Standalone chair (or linked to Eurus Cabinet Mounted Module or Cart) available in all seamless and luxury upholstery

£7,170

Eurus Contour PH Chair

with Power headrest (special order)
available in all seamless and luxury upholstery

£8,110

Please see swatch brochure for upholstery options

Folding & Fixed Legrest Dental Chairs



Eurus Swivel Chair

features include:

Standalone chair (or linked to Eurus Cabinet Mounted Module or Cart)
Maximum safe working patient weight of 200kg (ISO 7494-1)
Taller backrest cushion with additional lumbar support for greater patient comfort
Cable free base rotation 105 degrees (L & R)
Twin-articulating headrest with more precise positioning and lighter adjustment
Synchronised axis backrest, seat & legrest movement
Integrated foot stalks for chair control with 2 pre-set positions, Last Position Memory (LP) & Auto-return
Automatic extendable legrest with manual override & auto safety stop
Foldaway armrests
Seamless upholstery as standard

All standalone chairs can have a light post and bracket fitted (see page 43 for details)

Eurus Contour

features include:

Standalone chair (or linked to Eurus Cabinet Mounted Module or Cart)
Maximum safe working patient weight of 200kg (ISO 7494-1)
Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return
Swing-out armrest to 90 or 135 degrees
Synchronised axis backrest & seat movement
Integrated foot stalks for chair control
Seamless upholstery as standard

All Eurus Chairs are available with power headrest as a special order

Intra-Oral X-Ray Units

097 BELRAY II

097WK-300 Belray II wall-mounted, with 300mm primary arm, wall plate, scissor-type arm	£3,625
097WK-800 Belray II wall-mounted, with 800mm primary arm, wall plate, scissor-type arm	£3,625
097WK-1000 Belray II wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm	£3,625

Belmont Touch PHOT-X IIS DC LCD X-RAY

505WK-300 wall-mounted, with 300mm primary arm, wall plate, scissor-type arm	£4,305
505WK-800 wall-mounted, with 800mm primary arm, wall plate, scissor-type arm	£4,305
505WK-1000 wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm	£4,305
505FM mobile base	£4,740
505CK ceiling-mounted, with ceiling-mount pole & bracket	£5,195
505FK floor-mounted base	£4,030
505UM Unit-mounted for Eurus S4 Treatment Centre	£4,990

097 Belray II & Phot-X IIS DC LCD Intra-Oral X-Ray Units & optional extras



097 Belray II

features include:

70 KVP
0.7mm focal spot
Digital time display
Digital imaging compatible
(min. exposure 0.02 sec)
Manual override for total flexibility
Touch control for 2 film settings & digital imaging
Lockable touch pad control
Optional cone length & tube current
Rectangular Collimator
Hand exposure switch assembly

models include:

097WK-300

097WK-800

097WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach)
of 300mm (1,350*) / 800mm (1,850*) / 1,000mm
(2,050*)

Phot-X IIS DC LCD X-Ray

features include:

Touchscreen control
Easy-clean glass screen
Safety-conscious technologies
USB port
Instant Error Code reports
Lockable keypad
Adult/adolescent/child presets
60 or 70 KVP selectable
0.4mm focal spot
Digital time display
Radiation dose display
Compatible for digital imaging
(min. exposure 0.01 sec)
Manual override for total flexibility
Optional cone length, tube current & tube voltage
Rectangular Collimator
Hand exposure switch assembly

models include:

505WK-300

505WK-800

505WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach)
of 300mm (1,420*) / 800mm (1,920*) / 1,000mm
(2,120*)

505FM Mobile base

505FK Floor-mounted base

505CK Ceiling-mounted, with ceiling-mount
pole & bracket

* OPTIONS:

Ceiling Height (mm)	Pole Length (mm)
2,700	1,070
3,000	1,370

***Please specify when ordering**

Optional extras for 097 Belray II

£350	Long cone for paralleling technique
-------------	-------------------------------------

Optional extras for Phot-X IIS

£350	Long cone for paralleling technique
-------------	-------------------------------------

Operating Lights

100 LED SERIES MODELS

D-111 LED (Replaces 301U) post-mounted for Eurus standalone chairs only	£1,895
D-112 LED (Replaces 302C) ceiling-mounted	£2,740
D-119 LED (Replaces Voyager III 320LR) unit-mounted for Voyager III Plus only	£2,305

900 LED SERIES MODELS

D-901U LED post-mounted for Eurus standalone chairs	£2,965
D-902C LED ceiling-mounted, suitable for a maximum ceiling height 305cm	£3,295
D-902CX LED ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray	£3,395
D-905T LED ceiling-track - please specify pole length required	£4,290
D-921 LED unit-mounted, only compatible with Eurus S1, S4 & S6	£3,095
D-921LR LED unit-mounted, only compatible with Eurus S8	£3,095

REPLACEMENT LIGHTS FOR PREVIOUS UNITS

320LR LED + Kit unit-mounted light plus fitting kit 1P051PA0 (to replace the 048 & 720S LR lights for Voyager II/III only)	£2,445
920R LED + Kit unit-mounted light plus fitting kit for Cleo II 1P087VA0 unit-mounted light plus fitting kit for Clesta II 1P087VA0-CL2	£3,035
920LR LED + Kit unit-mounted light plus fitting kit for Voyager III 1P08CEA0 unit-mounted light plus fitting kit for tbCompass 1P08CDA0	£3,035
301U LED post-mounted light for 050, Progres and Clesta II chairs	£1,845

LED Operating Lights & optional extras



100 LED Series

features include:

The next generation 6 LED light
Reproduces natural light with a high functional colour rendering capacity
Central illuminance 5,000 to 34,000 lux
Correlated colour temperature 5000 K (4 LEDs)
Composite Curing Safe Mode (ISO9680 compliant)
Central illuminance: 11,000 lux
Correlated colour temperature: 2,700 K (2 LEDs)
Pattern dimensions 85 x 155mm
Colour rendering greater than 90
Irradiation distance 650mm
Touchless sensor switch for on/off and switching to composite mode
Auto ON/OFF with preset chair positions, unit mounted only
Triple-axis head positioning
One-piece removable front lens cover
Removable autoclavable operating handles
Optional patient viewing mirror

models include:

D-111 Post-mounted
D-112 Ceiling-mounted
D-119 Unit-mounted, only compatible with Voyager III Plus 'surgery system' & Voyager III Plus 'SP'

900 LED Series

features include:

The next generation 10 LED light
Reproduces natural light with a high functional colour rendering capacity
Central illuminance 5,500 to 39,000 lux
Correlated colour temperature 5000 K (8 LEDs)
Composite Curing Safe Mode (ISO9680 compliant)
Central illuminance: 7,500 lux
Correlated colour temperature: 2,800 K (2 LEDs)
Pattern dimensions 100x 170mm
Colour rendering greater than 90
Irradiation distance 650mm
Touchless sensor switch for on/off and switching to composite mode.
Auto ON/OFF with preset chair positions, unit mounted only
Triple-axis head positioning
Removable autoclavable operating handles
Optional patient viewing mirror

models include:

D-901U LED Post-mounted
D-902C LED Ceiling-mounted, suitable for a maximum ceiling height 305cm
D-902CX LED Ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray unit
D-905T LED Ceiling-track with pole:

* OPTIONS: Please specify when ordering

Ceiling Height (mm)	Pole Length (mm)
2,380	380
2,480	480
2,530	530
2,680	680
2,840	840
3,000	1,000

920LR LED Chair-mounted, only compatible with Voyager III and III 'SP'

D-921 LED Unit-mounted, only compatible with Eurus S1, S4 & S6

D-921LR LED Unit-mounted, only compatible with Eurus S8

Optional extras for Operating Lights

£915 Chair-mounted light-post and bracket for all chairs (Part No UK-217 & UK-218)

£215 1,500 mm ceiling-pole in lieu of 1,000 mm pole (suitable for max. ceiling height 3,550 mm) (Part No. UK-068)

£140 Patient Mirror for 100 & 900 LED Series Lights (Part No. D-900 1G02WMB0 D-100 1P05U9A0)

Dental Stools

SADDLE STOOL

Saddle Stool Gas Lift	£685
---------------------------------	-------------

FOR OPERATORS

ic-D Gas Lift	£475
ic-DF Gas Lift	£505
ic-Tilt Gas Lift	£665

FOR ASSISTANTS

ic-DSA Gas Lift	£515
ic-DSAF Gas Lift	£620
ic-DSSA Swing-Arm type	£715

Please see swatch brochure for upholstery options

Operators' & Assistants' Dental Stools & optional extras



Saddle Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Adjustable seat tilt angle
Seat height 580 - 780mm max.

ic-D Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Seat height 490 - 620mm max.

ic-DF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle
Seat height 500 - 630mm max.

ic-Tilt Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable seat & backrest angle
Seat height 500 - 630mm max.

ic-DSA Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Footrest ring
Seat height 530 - 730mm max.

ic-DSAF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle
Footrest ring
Seat height 540 - 740mm max.

ic-DSSA Stool

Swing-Arm type

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Backrest horizontal forward and back movement
Footrest ring
Seat height 590 - 730mm max.

* All stools benefit from:
Easy clean upholstery
Double-density cushioning

Optional extras for Operating Stools

£65 Footring for ic-Tilt, ic-DF & ic-D

Upholstery

Fabric range options for Treatment Centres, Dental Chairs and Dental Stools

SEAMLESS UPHOLSTERY

Standard colours	£0
Ultrasoft colours	£480
Ultrasoft Pearlised colours	£480

LUXURY UPHOLSTERY

Ultrasoft colours Not available on Voyager III treatment centre	£790
Ultrasoft Pearlised colours Not available on Voyager III treatment centres	£790

DENTAL STOOLS

Standard colours	£0
Ultrasoft colours	£105
Ultrasoft Pearlised colours	£105

IMPORTANT NOTE:

Cancelled or amended orders of Ultrasoft or Ultrasoft Pearlised upholstery will incur the following fee: £500 per Luxury chair, £250 per Seamless chair and £50 per stool.

Seamless Upholstery

Available in a choice of three fabrics



Eurus S1 Chair shown with Seamless Upholstery

Luxury Upholstery

Available in a choice of two fabrics



Eurus Swivel Chair shown with Luxury Upholstery

Eurus | S6



Our flagship folding legrest treatment centre allows **closer patient consultations** and **reduces patient anxiety** with discreet positioning of doctors table and instruments

 **Belmont**

Find out about our **new**
Dental Treatment Education Unit



The perfect start to a great dental career

See our brochure:



Hygiene & Care

For Cleaning Standards Beyond the Standard

For cleaning and maintenance of your treatment centre we highly recommend the following products:

TREATMENT CENTRE MAINTENANCE

B100 (250ml) Upholstery cleaning and care	£45
B300 (1 litre) Surface disinfectant	£38

Effective against all viruses, including COVID-19

Designed for easy use

Guarantees safety during application

Approved by the Commission for Hospital Hygiene & Infectious
Disease Prevention at the Robert Koch Institute (RKI)



Dental Equipment Service Kits

ANNUAL SERVICE KIT 1

£270

for : Voyager I, II, III & III Plus Treatment Centres

ANNUAL SERVICE KIT 2

£175

for : Voyager Module and Cart

ANNUAL SERVICE KIT 3

£190

for : Clesta I & II Module and Cart

ANNUAL SERVICE KIT 4

£305

for : Clesta II and Compass Treatment Centres

ANNUAL SERVICE KIT 5

£355

for : SP-Cleo and Cleo II Treatment Centres

ANNUAL SERVICE KIT 6

£355

for : Clesta I Treatment Centres

ANNUAL SERVICE KIT 7

£280

for : Eurus S1, S4, S6 & S8 Treatment Centres, Cabinet Mounted Module & Mobile Cart

Scan to download PDF price list:



SALES

LONDON & SOUTH EAST, THE SOUTH WEST, NORTHERN & SOUTHERN IRELAND

Eugene O'Malley **T:** 0771 097 9668

MIDLANDS, NORTHERN ENGLAND & SCOTLAND

Mark Harris **T:** 0786 053 5968

London Head Office & Showroom

Belmont House, One St. Andrews Way, Bow, London E3 3PA

T: 020 7515 0333

E: dental@takara.co.uk

Manchester Showroom

7 Quays Reach, Carolina Way, Salford, Manchester M50 2ZY

T: 0786 053 5968

E: dental@takara.co.uk

All showroom visits by appointment only

Belmont

belmontdental.co.uk