

A high-angle, close-up photograph of a dental chair. The chair is upholstered in a teal-colored material with visible stitching. To the left, a white dental handpiece is mounted on a control arm. Below the handpiece is a white dental sink. The background is a light, neutral color.

Ultrasoft Pro

colours for Luxury &
Seamless upholstery

Advanced fabric with antimicrobial
protection, ink and denim resistance

 **Belmont**

Ultrasoft Pro upholstery colours

for Belmont's Luxury and Seamless upholstery specifications

See product images in your preferred colours online at belmontdental.co.uk




Luxury Upholstery



Seamless Upholstery

SPECIFICATION

- Ink and denim resistant
- Performance plus
- Antimicrobial protection
- Custom
- Earth friendly



554-4315 Agave



554-0008 White Wash




554-6566 Moon Glow




554-5648 Cloud




554-5733 Mortar



554-2661 Blue Bird



554-2530 Sky



554-2522 Capri




554-4530 Coastal



554-4534 Sherwood



554-2711 Blueberry



554-2616 Oxford Blue



554-2547 Sailor



554-5626 Char



554-3811 Swan



554-3352 Shore



554-3595 Cork



554-5618 Dark Knight



554-3087 Woodcres



554-8251 Melon



554-8254 Rosy Coral



554-8229 Campfire



554-1329 Rhubarb



554-1357 Cranberry



554-9451 Azalea



554-9441 Mauve



Shown for representative purposes only.



Ultrasoft Pro

Ink and denim resistant, antimicrobial, soft-to-touch upholstery fabric

GENERAL SPECIFICATIONS

Composition:

100% Polycarbonate Polyurethane Surface
100% Reinforced Rayon Backing
containing TENCEL™ Modal fibers

Weight:

Approx. 9.7 oz. Per Square Yard
Approx. 330 Grams Per Square Metre

Thickness:

Approx. 0.95 Millimetres

Width:

54 Inches / 137 Centimetres

Repeat:

None

PERFORMANCE

Abrasion:

ASTM D4157 / 500,000 Double Rubs (Wyzenbeek)*
BS EN ISO 12947-2 / 200,000 Rubs (Martindale)*
*meets or exceeds

Hydrolysis:

14+ Weeks
ASTM D3690 Sect. 6.11 (Hydrolytic Stability)
ISO 1419 Method C (Accelerated Ageing)

Adhesion of Coating:	ASTM D751
Tear Strength:	ASTM D2261
Breaking Strength:	ASTM D751
Antimicrobial:	ASTM E2180
Tensile Strength:	BS EN ISO 1421

COLORFASTNESS

Light:	AATCC 16.3, 200 Hours, Grade 4 ISO 105 B02, Grade 6
---------------	--

Crocking:	AATCC 8, Grade 4-5 Wet & Dry
------------------	------------------------------

FLAMMABILITY REQUIREMENTS

- ASTM E84 (Adhered) - Class A
- CA TB 117-2013 - Pass
- FMVSS 302 & 95/28 Annex IV - Pass
- NFPA 260 - Class 1
- UFAC - Class 1

These fabrics, when treated and/or used with other compatible furniture components, have been shown to perform under the fire test conditions of:

- ASTM E84 (Unadhered)
- BS 5852 Source 0, 1 & 5 (Cig, Match & Crib 5)
- IMO FTPC Part 8
- BS EN 1021 Parts 1 & 2

Ultrasoft Pro Resists:

- | | |
|-----------------|---------------------|
| • Ballpoint ink | • Betadine |
| • Blood | • Denim* |
| • Red Wine | • Mustard & Ketchup |


* Removal of denim dye may vary depending on the period of time the stain sets in. While Ultrasoft Pro offers excellent protection, full removal is not guaranteed.

For more information on certifications and standards, sustainability, as well as cleaning and care instructions, please visit: ultrafabricsinc.com

Made in Japan



Perfecting the art of dentistry



For any enquiry, do not hesitate to contact us directly on 020 7515 0333 or visit: belmontdental.co.uk