

Dental Equipment **Retail Price List** April 2023 v2



belmontdental.co.uk

All prices listed are shown exclusive of VAT.

The manufacturers reserve the right to alter the prices and specifications without prior notice.

See our
website:



Tell me about the
Eurus touch-screen:



CONTENTS

- 2 **FREE 5 Year Warranty**
- 4 **Eurus S1**
Treatment Centres / Flexible Options / Factory Built-in Options
- 10 **Eurus S6**
Treatment Centres / Flexible Options / Factory Built-in Options
- 16 **Eurus S8**
Treatment Centres / Factory Built-in Options
- 22 **Voyager III**
Treatment Centres / Flexible Options / Factory Built-in Options
- 26 **Flexible Delivery Options**
Cleo II Flexible / Voyager 'M' / Voyager 'MC' / Clesta II 'M' / Clesta II 'MC'
- 30 **Dental Chairs**
Progres / Clesta II / 050 Clair / 055 Clair (Special Order) /
Eurus Swivel / Eurus Swivel PH (Special Order)
- 32 **Intra Oral X-rays**
097 Belray II / Belmont Touch Phot-X IIS DC LCD
- 34 **Dental Operating Lights**
300 LED Series / 900 LED Series
- 36 **Dental Stools**
For Operators & Assistants
- 38 **Upholstery**
Standard / Ultrasoft / Ultrasoft Pro / Pearlised
- 40 **Hygiene & Care**
Cleaning and maintenance range

FREE

5 Year Warranty



Takara Belmont (UK) Ltd Warranty Policy

5 Year Warranted Products

All Belmont treatment centres, individual chairs, units and dental lights are warranted for a period of five years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions.

Upholstery, paintwork, plating and all handpiece hoses are not warranted, unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery. Failure due to natural wear and tear of any component is not warranted.

OEM factory build-in options (ie: fibre optics, scalers, motors, separators, CVS Kits etc) carry a warranty as recommended by the individual manufacturer.

2 Year Warranted Products

All Belmont intra oral x-rays are warranted for two years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions. Paintwork and plating are not warranted unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery.

Takara Belmont (UK) Ltd are not responsible for any failure caused by misuse of the equipment or damage caused due to poor installation by a non-authorised company.

The warranty policy applies to equipment purchased from and installed by an authorised Takara Belmont (UK) Ltd dealer. Labour will be charged on warranty work after the first year.

Hygiene & Care

In order to maintain & prolong the longevity of your Belmont equipment, we highly recommend using our Hygiene & Care products. (See page 40)

PLEASE NOTE: ISO 13485 Certification

All of our dental products are supplied with a full CE mark and certification to the standard complying with ISO 13485. Any modifications or alterations made after certification may affect your warranty as the responsibility for CE and electrical safety will transfer to the modifier.

TREATMENT CENTRES (Over-the-Patient delivery system)	With Unit Mounted 921 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S1 'Surgery System' Holder Type (A)	£20,970	£21,410	£19,320
Eurus S1 'Surgery System' Holder Type (E)	£23,830	£24,270	£22,180
Eurus S1 'Surgery System' CM Rod Type (A)	£22,160	£22,600	£20,510
Eurus S1 'Surgery System' CM Rod Type (E)	£25,020	£25,460	£23,370
Eurus S1 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£21,495	£21,935	£19,845
Eurus S1 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£24,355	£24,795	£22,705
When light sold with package above the saving over recommended retail price is:	Save £1,190	Save £935	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Chair Mounted Over-Patient Holder & Rod Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation

consists of ~

Chair section features include:

The Eurus S1 has a maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Integrated service centre
Synchronised axis backrest & seat movement
Integrated foot stalks for chair control
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Autoclavable unit handle cover, handpiece holders
unit tray silicon mat, handpiece rest (Rod Type)
Easy-clean membrane switch panel
Doctor handpiece control touch screen with water low warning indicator
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Electronic foot control for accurate turbine speed control and integrated lifting handle

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system
2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation
921 LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Bottled water system with water low warning
LED on cup filler spout casing
Durr spittoon valve



Eurur S1 Flexible System
with Mobile Cart

Eurur S1 Treatment Centres



Flexible Treatment Centres Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation
consists of ~

Chair section features include:

The Eurur S1 has a maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Integrated service centre
Synchronised axis backrest & seat movement
Integrated foot stalks for chair control
Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurur Cabinet Module or Mobile Cart delivery systems A & E (please specify)

'A' version features include:

Cabinet mounting bracket and arm or height adjustable mobile cart base
Autoclavable unit handle cover, handpiece holders and unit tray silicon mat
Easy clean membrane switch panel
Doctor handpiece control touch screen
3-way autoclavable syringe
3 Fibre optic & 1 non optic handpiece outlets with individual air & water controls
First priority handpiece selection
Electronic foot control for accurate turbine speed control and integrated lifting handle
Handpiece flush-out
Handpiece oil mist collector
Bottled water system

'E' version additional features include:

Bien-Air MX2 fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system
2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation
921 LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Bottled water system with water low warning
LED on cup filler spout casing
Durr spittoon valve

Important:

Maximum loom length between module or cart to the chair is 5.7metres (standard is 2.7metres).

Eurus S1



Eurus S1 shown with Luxury Upholstery

Eurus S1 Factory Built-in Options



for Eurus S1	
£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,375	Additional Bien Air MX2 fibre-optic micro motor *
£1,640	Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence
£1,535	Additional NSK NLX-Plus fibre-optic micro motor *
£1,995	Cavitron G-139 scaler
£1,405	NSK Varios 170 non-optic scaler
£1,670	NSK Varios 170 optic scaler
£1,470	Satelec Newtron scaler
£1,735	Satelec Newtron LED scaler
£980	Satelec mini LED curing light
£2,775	Durr CAS 1 Combi-Separator with rinse
£3,690	Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount)
£4,935	Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount)
£905	Wireless unit foot controller (when factory fitted)
£680	Monitor mount (if not supplied with camera kit)
£390	Additional 3-way autoclavable syringe (assistant side)
£190	Additional easily detachable HV ejector hose (including ball joint for greater ease of operation)
£175	Large tray (double size for Rod Type)
£85	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£225	Flushing Tray A (Holder Type)
£225	Flushing Tray B (Rod Type)

* Additional micro motors must be the same model as the primary micro motor.



New Flushing Tray - suitable for all treatment centres

TREATMENT CENTRES (Below-the-Patient delivery system)	With Unit Mounted 921 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S6 'Surgery System' Holder Type (A)	£26,550	£26,990	£24,900
Eurus S6 'Surgery System' Holder Type (E)	£29,410	£29,850	£27,760
Eurus S6 'Vac-Pac System' Holder Type (A)	£22,030	£22,470	£20,380
Eurus S6 'Vac-Pac System' Holder Type (E)	£24,890	£25,330	£23,240
Eurus S6 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£25,880	£26,320	£24,230
Eurus S6 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System	£28,740	£29,180	£27,090
When light sold with package above the saving over recommended retail price is:	Save £1,190	Save £935	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Eurus S6 Treatment Centres with folding and extending legrest



Surgery System (A and E)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Taller backrest cushion with additional lumbar support for greater patient comfort

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray with pneumatic release and lock button

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders unit tray silicon mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate

instrument operation with integrated lifting handle

Bottle water system

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system 2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with

90 degree rotation

Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left side 921 LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the cup filler spout and touch-screen

Durr spittoon valve

'Vac-Pac System' (A and E)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

Assistant's console

Easy-access detachable filter

Autoclavable HVE & saliva ejector

Funnel connection hose



Eurus S6 enables greater workflow freedom

Eurus S6 Treatment Centres with folding and extending legrest



Flexible Treatment Centres Type (A and E)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg
(IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise
positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) &
auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the
back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Taller backrest cushion with additional lumbar
support for greater patient comfort

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurus Cabinet Module or
Mobile Cart delivery systems A & E (please specify)

'A' version features include:

Cabinet mounting bracket and arm or height
adjustable mobile cart base

Autoclavable unit handle cover, handpiece holders
and unit tray silicon mat

Easy clean membrane switch panel

Doctor handpiece control touch screen

3-way autoclavable syringe

3 Fibre optic & 1 non optic handpiece outlets with
individual air & water controls

First priority handpiece selection

Wireless electronic foot control for accurate

instrument operation with integrated lifting handle

Handpiece flush-out

Handpiece oil mist collector

Bottled water system

'E' version additional features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with

auto-reverse, Endo control and LCD display, with
optional up-grade to WaveOne endodontic system

2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with
90 degree rotation

Cuspidor with 90 degree rotation away from the
chair for accessing the patient from the left side

921 LED dental light with detachable autoclavable
handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the
cup filler spout and touch-screen

Durr spittoon valve

Important:

**Maximum loom length between module
or cart to the chair is 5.7metres
(standard is 2.7metres).**

Eurus S6



Eurus S6 shown with Luxury Upholstery

Eurus S6 Factory Built-in Options



for Eurus S6	
£0	NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only)
£1,375	Additional Bien Air MX2 fibre-optic micro motor *
£1,640	Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence
£1,535	Additional NSK NLX-Plus fibre-optic micro motor *
£1,995	Cavitron G-139 scaler
£1,405	NSK Varios 170 non-optic scaler
£1,670	NSK Varios 170 optic scaler
£1,470	Satelec Newtron scaler
£1,735	Satelec Newtron LED scaler
£980	Satelec mini LED curing light
£2,775	Durr CAS 1 Combi-Separator with rinse
£3,690	Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount)
£4,935	Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount)
£680	Monitor mount (if not supplied with camera kit)
£390	Additional 3-way autoclavable syringe (assistant side)
£85	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£225	Flushing Tray A (Holder Type)

* Additional micro motors must be the same model as the primary micro motor.

TREATMENT CENTRES (Below-the-Patient delivery system)	With Unit Mounted 921LR Light	With Ceiling Mounted D-902C Light	Without Dental Light
Eurus S8 'Surgery System' Holder Type (A)	£22,920	£23,360	£21,270
Eurus S8 'Surgery System' Holder Type (E)	£25,780	£26,220	£24,130
Eurus S8 'Vac-Pac System' Holder Type (A)	£20,730	£21,170	£19,080
Eurus S8 'Vac-Pac System' Holder Type (E)	£23,590	£24,030	£21,940
Eurus S8 'Surgery System' CM Rod Swing (A)	£25,935	£26,375	£24,285
Eurus S8 'Surgery System' CM Rod Swing (E)	£28,795	£29,235	£27,145
Eurus S8 'Vac-Pac System' CM Rod Swing (A)	£23,065	£23,505	£21,415
Eurus S8 'Vac-Pac System' CM Rod Swing (E)	£25,925	£26,365	£24,275
When light sold with package above the saving over recommended retail price is:	Save £1,190	Save £935	£0

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg
(IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Synchronised axis backrest & seat movement
Integrated foot stalk or combined with unit foot controller (Holder Type)
Integrated foot stalks for chair control (Rod Type)
Integrated service centre
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system
Height adjustable unit tray with 3 settings (Holder Type)
Powerful intuitive touch-screen with advanced instrument control
Autoclavable unit handle cover, handpiece holders
unit tray silicon mat
Easy-clean membrane switch panel
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Wireless electronic foot control for accurate instrument operation with integrated lifting handle plus chair manual control movements (Holder Type)
Bottle water system

Unit section (E) Doctor side additional

features include:

NSK NLX-Plus fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse & touch-screen display
2 fibre-optic outlets

WaveOne is available as an additional option when the Bien Air MX2 micro motor is fitted

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 90 degree rotation
Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left and right sides
921LR LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
2 x HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Chair movement lock indicators located on the cup filler spout and touch-screen
Durr spittoon valve



Base mounted ambidextrous
Holder-Type (A and E) treatment centre

Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg
(IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Synchronised axis backrest & seat movement
Integrated foot stalk combined with unit foot controller (Holder Type)
Integrated foot stalks for chair control (Rod Type)
Integrated service centre
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system
Height adjustable unit tray with 3 settings (Holder Type)
Powerful intuitive touch-screen with advanced instrument control, plus WaveOne technology option
Autoclavable unit handle cover, handpiece holders & tray silicon mat / Rod Type handpiece place mat
Easy-clean membrane switch panel
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Wireless electronic foot control for accurate instrument operation with integrated lifting handle plus chair manual control movements (Holder Type)
Bottle water system

Unit section (E) Doctor side additional

features include:

NSK NLX-Plus fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse & touch-screen display
2 fibre-optic outlets

'Vac-Pac System' (A and E)

(Ambidextrous)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-access detachable filter
Autoclavable HVE & saliva ejector
Funnel connection hose
Chair control auto & manual touch pad
Assistant tray

Eurus S8

Chair mounted under patient ambidextrous
Rod Swing (A and E) Treatment Centre shown
with optional large tray



Eurus S8 shown with Luxury Upholstery

Eurus S8 Factory Built-in Options



for Eurus S8	
£0	Bien Air MX2 micro motor (in place of NSK NLX-Plus micro motor on 'E' version only)
£1,375	Additional Bien Air MX2 fibre-optic micro motor *
£1,640	Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence
£1,535	Additional NSK NLX-Plus fibre-optic micro motor *
£1,995	Cavitron G-139 scaler
£1,405	NSK Varios 170 non-optic scaler
£1,670	NSK Varios 170 optic scaler
£1,470	Satelec Newtron scaler
£1,735	Satelec Newtron LED scaler
£980	Satelec mini LED curing light
£2,775	Durr CAS 1 Combi-Separator with rinse
£3,690	Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount)
£4,935	Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount)
£680	Monitor mount (if not supplied with camera kit)
£390	Additional 3-way autoclavable syringe (assistant side)
£175	Large tray (double size for Rod Type)
£85	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£225	Flushing Tray A (Holder Type)
£225	Flushing Tray B (Rod Type)

* Additional micro motors must be the same model as the primary micro motor.

TREATMENT CENTRES Below-the-Patient delivery system and flexible options	With Unit Mounted 320LR Light	With Ceiling Mounted 302C Light	Without Dental Light
Voyager III 'Surgery System'	£14,635	£15,020	£13,350
Voyager III 'SP' Chair with ambidextrous cuspidor with Voyager cabinet mounted module or mobile cart delivery system	£14,635	£15,020	£13,350
When light sold with package above the saving over recommended retail price is:	Save £825	Save £835	£0

TREATMENT CENTRES Below-the-Patient delivery system and flexible options	With Unit Mounted 920LR Light	With Ceiling Mounted D-902C Light	Without Dental Light
Voyager III 'Surgery System'	£15,190	£15,970	£13,350
Voyager III 'SP' As above	£15,190	£15,970	£13,350
When light sold with package above the saving over recommended retail price is:	Save £825	Save £935	£0

Voyager III 'SP' Chair with ambidextrous cuspidor (without operating light)	£10,210
--	---------

Please state light model number when ordering treatment centre

MODULE & MOBILE CART OPTIONS

Voyager 'M' & 'MC' Module and Mobile Cart delivery system	£3,140
--	--------

Voyager III Treatment Centres & Flexible Systems

'Surgery System'

(Ambidextrous)

with 320LR unit-mounted Operating Light

consists of ~

Chair section features include:

The Voyager III has a maximum safe working patient weight of 135kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrests (left & right)
Integral service centre with air & water filter regulators
8-way foot control

Base-mounted swing-around delivery system features include:

3 handpiece outlets with individual air & water control
3-way autoclavable syringe
Table safety lock device
Quick flush-out system
Bottled water system
Stainless steel instrument tray
Non-slip mat

Chair-mounted rotatable cuspidor unit

features include:

Porcelain cuspidor bowl
Manual cup-filler & bowl-flush
1 detachable high volume evacuator
1 detachable saliva ejector
Easily-accessed & detached solids collector
Durr spittoon valve

'SP'

(Ambidextrous)

Voyager III Dental Chair with cuspidor

consists of ~

Chair section features include:

The Voyager III has a maximum safe working patient weight of 135kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrests (left & right)
Integral service centre with air & water filter regulators
8-way foot control

Chair-mounted rotatable cuspidor unit features include:


Porcelain cuspidor bowl
Manual cup-filler & bowl-flush
1 detachable high volume evacuator
1 detachable saliva ejector
Easily-accessed & detached solids collector
Durr spittoon valve

Voyager III Factory Built-in Options

for Voyager

£2,405	Bien Air Optima MCX LED (inc. standard fibre-optic micro motor, turbine hose & touch pad control)
£2,405	NSK NLX-Nano fibre-optic pack (inc. NLX-Nano fibre-optic micro motor, turbine hose & touch pad control)
£2,935	NSK NLX-Plus multi-pad endo pack (inc. NLX-Plus fibre-optic micro motor, turbine hose & multi-pad Endo control)
£835	Bien Air fibre-optic kit (inc. 1 turbine hose)
£285	Additional fibre-optic outlet
£2,105	Cavitron G-139 scaler
£1,405	NSK Varios 170 non-optic basic scaler
£1,670	NSK Varios 170 optic scaler (with touch pad control)
£1,505	Satelec Newtron scaler
£1,780	Satelec Newtron LED scaler
£1,170	Satelec mini LED curing light
£420	Additional 3-way autoclavable syringe on Assistant's side inc. bottled water system
£0	Standard (10x10) Module Arm (alternative to Voyager Arm System)
£495	Touch pad assembly for chair control fitted to Operator's table
£795	Monitor mount (if not supplied with camera kit)
£3,690	Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount)
£4,935	Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount)
£80	Extra-long umbilical (up to 4 metres for Mobile Cart only)
£85	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£115	Floor box fitted to cart base
£225	Flushing Tray A (Holder Type)

Did you know you can have a private viewing?



**Our showrooms are open during office hours,
weekends & evenings, whilst our 3D showroom
experience is available from the comfort of your home.**

You can access our showroom 24/7 from wherever you are,
simply go to belmontdental.co.uk/showroom

Flexible Delivery Options

CLEO II FLEXIBLE TREATMENT CENTRES	With Unit Mounted 920 Light	With Ceiling Mounted D-902C Light	Without Dental Light
Cleo II 'Flexible System' (A) Dental chair and cuspidor with Clesta II cabinet mounted module or mobile cart delivery system	£22,945	£23,610	£21,520
Cleo II 'Flexible System' (E) Dental chair and cuspidor with Clesta II cabinet mounted module or mobile cart delivery system	£26,285	£26,950	£24,860
When light sold with package above the saving over recommended retail price is:	Save £1,240	Save £935	£0

CHAIR & CUSPIDOR OPTION (without Operating Light)

Cleo II 'Flexible System' dental chair and cuspidor	£16,525
---	---------

MODULE & MOBILE CART OPTIONS

Voyager 'M' & 'MC' Module and Mobile Cart delivery system	£3,140
Clesta II 'M' & 'MC' (A) Module and Mobile Cart delivery system	£4,995
Clesta II 'M' & 'MC' (E) Module and Mobile Cart delivery system	£7,895
Clesta II 'M' (E) or Clesta II 'MC' (E) Module or Mobile Cart delivery system when supplied with Cleo II 'Flexible' treatment centre (E)	£8,335

Clesta II 'M' & 'MC' **Factory Built-in Options**
For Built-in Options please refer to Clesta prices on page 29

Voyager 'M' & 'MC' **Factory Built-in Options**
Please refer to page 24 for prices

Universally compatible Flexible Options



'Flexible System'

Cleo II Dental Chair and cuspidor consists of ~

Chair section features include:

Twin-articulating headrest
2 preset positions, Last-Position memory (LP) & Auto-Return
Automatic extendable legrest with manual override & Auto Safety Stop
Foldaway, detachable armrest
Integral service centre with air & water regulators & pressure indicators

Chair-mounted cuspidor features include:

Rotatable cuspidor (90°)
Assistant's console with balanced positioning
Assistant 3-way autoclavable syringe
Easily-detached HV & saliva ejector tubings
Automatic bowl-flush & cup-filler
Bottled water system
Assistant-side chair preset controls
Porcelain cuspidor bowl
Durr spittoon valve
~ **Note:** cuspidor functionality will require upgrading at additional cost if a Clesta II 'M' or 'MC' (E) is added post-delivery

Voyager 'M'

Module delivery system features include:

Mounting arms & bracket
3 handpiece outlets with individual air & water control
3-way autoclavable syringe
Bottled water system
Stainless steel instrument tray
Non-slip mat
Service centre with air filter regulator

Voyager 'MC'

Mobile Cart delivery system includes the same as Voyager 'M' but with height-adjustable mobile base

Clesta II 'M' (A)

Module delivery system features include:

Mounting arms & bracket
3-way autoclavable syringe
3 handpiece outlets with individual air & water control
Handpiece flush-out system
Bottled water system
Service centre with air filter regulator

Clesta II 'MC' (A)

Mobile Cart delivery system includes the same as Clesta II 'M' (A) but with height-adjustable mobile base

Clesta II 'M' (E)

Module delivery system includes the same as the 'A' type but with the following additional features:

NSK NLX-Plus fibre-optic micro motor
Preset micro motor & digital display
2 fibre-optic outlets
First Priority handpiece selection
Handpiece oil mist collector
Additional chair preset positioning

Clesta II 'MC' (E)

Mobile Cart delivery system includes the same as Clesta II 'M' (E) but with height-adjustable mobile base

Individual Delivery Options – p.26 Intra-Oral X-Ray Units – p.32 Operating Lights – p.34
Upholstery Prices – p.38 Upholstery Options – see separate swatch brochure for colours



for Cleo II & Clesta II

£0	Bien Air MX2 micro motor (in place of NSK NLX-Plus micro motor on 'E' version only)
£3,475	Bien Air MX2 Optima INT Endo pack (inc. standard fibre-optic micro motor, turbine hose & Optima INT Endo control)
£2,380	NSK NLX-Nano fibre-optic pack (inc. NLX-Nano fibre-optic micro motor & turbine hose)
£2,935	NSK NLX-Plus multi-pad Endo pack (inc. NLX-Plus fibre-optic micro motor, turbine hose & multi-pad Endo control)
£2,045	Additional Bien Air MX2 fibre-optic micro motor *
£1,075	Additional NSK NLX-Nano fibre-optic micro motor *
£1,535	Additional NSK NLX-Plus fibre-optic micro motor *
£775	Bien Air fibre-optic kit (inc. 1 turbine hose)
£285	Additional fibre-optic outlet
£1,995	Cavitron G-139 scaler
£1,405	NSK Varios 170 non-optic scaler
£1,670	NSK Varios 170 optic scaler
£1,470	Satelec Newtron scaler
£1,735	Satelec Newtron LED scaler
£980	Satelec mini LED curing light
£3,690	Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount)
£4,935	Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount)
£160	Handpiece oil mist collector (Flexible 'A' version only)
£80	Extra-long umbilical (up to 4 metres for Mobile Cart only)
£115	Floor box fitted to cart base
£2,690	Durr CAS 1 Combi-Separator without rinse
£2,775	Durr CAS 1 Combi-Separator with rinse
£795	Monitor mount (if not supplied with camera kit)
£85	Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm)
£225	Flushing Tray A (Holder Type)

Dental Chairs

WITH FOLDING AND EXTENDING LEGREST

Eurus Swivel Chair Standalone chair with Twin-articulating headrest & cable free base rotation 105 degrees (L & R) available in all seamless and luxury upholstery	£9,455
Eurus Swivel PH Chair with Power headrest (special order) available in all seamless and luxury upholstery	£10,485
050 Clair Chair with Twin-articulating headrest available in all seamless and luxury upholstery	£8,995
055 Clair Chair with Power headrest (special order) available in all seamless and luxury upholstery	£9,945

WITH FIXED LEGREST

Clesta II Chair Important note: to be used as standalone only	£6,760
Progres Chair features 60° (L & R) base rotation Not to be used with an Over-Patient unit Important note: available with Seamless upholstery only	£5,895

Please see swatch brochure for upholstery options

Folding & Fixed Legrest Dental Chairs



Eurus Swivel Chair

features include:

Standalone chair
Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)
Cable free base rotation 105 degrees (L & R)
Twin-articulating headrest with more precise positioning and lighter adjustment
Synchronised axis backrest, seat & leg-rest movement
Integrated foot stalks for chair control with 2 pre-set positions, Last Position Memory (LP) & Auto-return
Automatic extendable leg-rest with manual override & auto safety stop
Foldaway armrests
Seamless upholstery as standard

Eurus Swivel PH Chair (special order)

features include the same as above but with a Power headrest

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

050 Clair Chair

features include:

Maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Automatic extendable legrest with manual override & auto-safety stop
Foldaway, detachable armrests
Seamless upholstery

055 Clair Chair (special order)

features include the same as 050 Clair Chair but with a Power headrest

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

Clesta II Chair

features include:

Maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrest
Synchronised axis backrest
8-way foot control
Seamless upholstery

Important note: to be used as standalone only

Progres Chair

features include:

Maximum safe working patient weight of 135kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
Pantograph base with 60° rotation
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrest
8-way foot control
Seamless upholstery
The seating area can be rotated 30 degrees to the left or right and then locked into position

Important note: not to be used with an Over-Patient unit

Intra-Oral X-Ray Units

097 BELRAY II

097WK-300 Belray II wall-mounted, with 300mm primary arm, wall plate, scissor-type arm	£3,330
097WK-800 Belray II wall-mounted, with 800mm primary arm, wall plate, scissor-type arm	£3,330
097WK-1000 Belray II wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm	£3,330

Belmont Touch PHOT-X IIS DC LCD X-RAY

505WK-300 wall-mounted, with 300mm primary arm, wall plate, scissor-type arm	£3,965
505WK-800 wall-mounted, with 800mm primary arm, wall plate, scissor-type arm	£3,965
505WK-1000 wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm	£3,965
505FM mobile base	£4,360
505CK ceiling-mounted, with ceiling-mount pole & bracket	£4,780
505FK floor-mounted base	£3,705

097 Belray II & Phot-X IIS DC LCD Intra-Oral X-Ray Units & optional extras



097 Belray II

features include:

70 KVP
0.7mm focal spot
Digital time display
Digital imaging compatible
(min. exposure 0.02 sec)
Manual override for total flexibility
Touch control for 2 film settings & digital imaging
Lockable touch pad control
Optional cone length & tube current
Rectangular Collimator
Hand exposure switch assembly

models include:

097WK-300

097WK-800

097WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach) of
300mm (1,350*) / 800mm (1,850*) / 1,000mm
(2,050*)

Phot-X IIS DC LCD X-Ray

features include:

Touchscreen control
Easy-clean glass screen
Safety-conscious technologies
USB port
Instant Error Code reports
Lockable keypad
Adult/adolescent/child presets
60 or 70 KVP options
0.4mm focal spot
Digital time display
Radiation dose display
Compatible for digital imaging
(min. exposure 0.01 sec)
Manual override for total flexibility
Optional cone length, tube current & tube voltage
Rectangular Collimator
Hand exposure switch assembly

models include:

505WK-300

505WK-800

505WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach) of
300mm (1,420*) / 800mm (1,920*) / 1,000mm
(2,120*)

505FM Mobile base

505FK Floor-mounted base

505CK Ceiling-mounted, with ceiling-mount
pole & bracket

*** OPTIONS:**

Ceiling Height (mm)	Pole Length (mm)
2,700	1,070
3,000	1,370
3,200	1,570

***Please specify when ordering**

Optional extras for 097 Belray II

£310	Long cone for paralleling technique
-------------	-------------------------------------

Optional extras for Phot-X IIS

£310	Long cone for paralleling technique
£330	Operating Light bracket for 505CK

Operating Lights

300 LED SERIES MODELS

301U LED post-mounted	£1,735
302C LED ceiling-mounted	£2,505
320LR LED unit-mounted	£2,110
320LR LED + Kit unit-mounted light plus fitting kit 1P051PA0-KIT (to replace the 048 & 720S LR lights for Voyager II/III only)	£2,110

900 LED SERIES MODELS

901 LED post-mounted	£2,725
D-902C LED ceiling-mounted, suitable for a maximum ceiling height 305cm	£3,025
D-902CX LED ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray	£2,780
905 LED ceiling-track - please specify pole length required	£3,940
920 LED unit-mounted, only compatible with Cleo II	£2,665
920LR LED chair-mounted, only compatible with Voyager III & III 'SP'	£2,665
921 LED unit-mounted, only compatible with Eurus S1 & S6	£2,840
921LR LED unit-mounted, only compatible with Eurus S8	£2,840

LED Operating Lights & optional extras



300 LED Series

features include:

White natural beam, 5000 Kelvin (5 LEDs)
Touchless sensor switch
Composite Curing Safe Mode
Triple-axis head movement
Auto ON/OFF with preset chair positions*

models include:

301U LED Post-mounted
302C LED Ceiling-mounted
320LR LED Unit-mounted, *only compatible with Voyager III 'Surgery System' & Voyager III 'SP'*

* Applies to unit-mounted versions only.

900 LED Series

features include:

White natural beam, 5000 Kelvin (10 LEDs)
Touchless sensor switch
Composite Curing Safe Mode
Triple-axis head movement
Auto ON/OFF with preset chair positions

models include:

901 LED Post-mounted
D-902C LED Ceiling-mounted, suitable for a maximum ceiling height 305cm
D-902CX LED Ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray unit
905 LED Ceiling-track with pole:

* OPTIONS:

Ceiling Height (mm)	Pole Length (mm)
2,380	380
2,480	480
2,530	530
2,680	680
2,840	840
3,000	1,000

920LR LED Chair-mounted, only compatible with Voyager III and III 'SP'

921 LED Unit-mounted, only compatible with Eurus S1 & S6

921LR LED Unit-mounted, only compatible with Eurus S8

***Please specify when ordering**

Optional extras for Belmont Operating Lights

£535	Chair-mounted light-post & bracket for 055 & 050 Clair standalone chairs (Part No. UK-149)
£535	Chair-mounted light-post & bracket for Progres chairs (Part No. UK-194)
£195	1,500 mm ceiling-pole in lieu of 1,000 mm pole (suitable for max. ceiling height 3,550 mm) (Part No. UK-068)
£125	Patient Mirror for 900 LED Series (Part No. 1G02WMB0)

Dental Stools

SADDLE STOOL

Saddle Stool Gas Lift	£595
---------------------------------	-------------

FOR OPERATORS

ic-D Gas Lift	£415
ic-DF Gas Lift	£440
ic-Tilt Gas Lift	£580

FOR ASSISTANTS

ic-DSA Gas Lift	£450
ic-DSAF Gas Lift	£540
ic-DSSA Swing-Arm type	£625

Please see swatch brochure for upholstery options

Operators' & Assistants' Dental Stools & optional extras



Saddle Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Adjustable seat tilt angle

Seat height 580 - 780mm max.

ic-D Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height

Seat height 490 - 620mm max.

ic-DF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle

Seat height 500 - 630mm max.

ic-Tilt Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable seat & backrest angle

Seat height 500 - 630mm max.

ic-DSA Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Footrest ring

Seat height 530 - 730mm max.

ic-DSAF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle
Footrest ring

Seat height 540 - 740mm max.

ic-DSSA Stool

Swing-Arm type

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Backrest horizontal forward and back movement
Footrest ring

Seat height 590 - 730mm max.

* All stools benefit from:
Easy clean upholstery
Double-density cushioning

Optional extras for Operating Stools

£50	Footring for ic-Tilt , ic-DF & ic-D
------------	-------------------------------------

Upholstery

Fabric range options for Treatment Centres, Dental Chairs and Dental Stools

SEAMLESS UPHOLSTERY

Standard colours	£0
Ultrasoft colours	£395
Ultrasoft Pro colours	£395
Ultrasoft Pearlised colours	£395

LUXURY UPHOLSTERY

Ultrasoft colours Not available on Voyager III treatment centres & Progres chairs	£675
Ultrasoft Pro colours Not available on Voyager III treatment centres & Progres chairs	£675
Ultrasoft Pearlised colours Not available on Voyager III treatment centres & Progres chairs	£675

DENTAL STOOLS

Standard colours	£0
Ultrasoft colours	£90
Ultrasoft Pro colours	£90
Ultrasoft Pearlised colours	£90

IMPORTANT NOTE:

Cancelled or amended orders of Ultrasoft, Ultrasoft Pro or Ultrasoft Pearlised upholstery will incur the following fee: £500 per Luxury chair, £250 per Seamless chair, and £50 per stool.

Seamless Upholstery

Available in a choice of four fabrics



Eurus S1 Chair shown with Seamless Upholstery

Luxury Upholstery

Available in a choice of three fabrics



Eurus Swivel Chair shown with Luxury Upholstery

Hygiene & Care

For Cleaning Standards Beyond the Standard

For cleaning and maintenance of your treatment centre we highly recommend the following products:

TREATMENT CENTRE MAINTENANCE

B100 (500ml) Upholstery cleaning and care	£35
B300 (1 litre) Surface disinfectant	£30
B500 DUWL Starter Kit - 2 periodic shock tests, 2 periodic shock treatments and 20 daily maintenance sachets	£52

DENTAL UNIT WATERLINE PROTECTION

B700 (100 sachets) DUWL daily maintenance	£95
B900 (5 x 200ml) DUWL periodic shock and test treatment	£75

Effective against all viruses, including COVID-19

Designed for easy use

Guarantees safety during application

Approved by the Commission for Hospital Hygiene & Infectious Disease Prevention at the Robert Koch Institute (RKI)

We're offering a **free Intro Kit** with every new treatment centre ordered. This includes a **B100**, **B300** and **B500** as an introduction to our hygiene & care products.



Scan to download PDF price list:



SALES

LONDON & SOUTH EAST

Eugene O'Malley **T:** 0771 097 9668

MIDLANDS, SOUTH WEST, NORTHERN & SOUTHERN IRELAND

Barry Hunt **T:** 0771 526 3587

NORTHERN ENGLAND & SCOTLAND

Mark Harris **T:** 0786 053 5968

London Head Office & Showroom

Belmont House, One St. Andrews Way, Bow, London E3 3PA

T: 020 7515 0333

E: dental@takara.co.uk

Manchester Showroom

7 Quays Reach, Carolina Way, Salford, Manchester M50 2ZY

T: 0161 743 9992

E: dental@takara.co.uk

All showroom visits by appointment only

Belmont

belmontdental.co.uk