

Dental Equipment Retail Price List

April 2024



belmontdental.co.uk

All prices listed are shown exclusive of VAT.

The manufacturers reserve the right to alter the prices and specifications without prior notice.

See our
website:



Tell me about the
Eurus touch-screen:



CONTENTS

- 2 **FREE 5 Year Warranty**
- 4 **Eurus S1**
Treatment Centres / Flexible Options / Factory Built-in Options
- 10 **Eurus S6**
Treatment Centres / Flexible Options / Factory Built-in Options
- 16 **Eurus S8**
Treatment Centres / Factory Built-in Options
- 22 **Voyager III**
Treatment Centres / Flexible Options / Factory Built-in Options
- 26 **Flexible Delivery Options**
Voyager 'M' / Voyager 'MC' / Clesta II 'M' / Clesta II 'MC'
- 30 **Dental Chairs**
Clesta II / 050 Clair / 055 Clair /
Eurus Swivel / Eurus Swivel PH (Special Order)
- 32 **Intra Oral X-rays**
097 Belray II / Belmont Touch Phot-X IIS DC LCD
- 34 **Dental Operating Lights**
300 LED Series / 900 LED Series
- 36 **Dental Stools**
For Operators & Assistants
- 38 **Upholstery**
Standard / Ultrasoft / Ultrasoft Pro / Pearlised
- 41 **Hygiene & Care**
Cleaning and maintenance range
- 42 **Service Kits**
For Treatment Centres, Modules & Mobile Carts

FREE 5 Year Warranty



Takara Belmont (UK) Ltd Warranty Policy

5 Year Warranted Products

All Belmont treatment centres, individual chairs, units and dental lights are warranted for a period of five years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions.

Upholstery, paintwork, plating and all handpiece hoses are not warranted, unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery. Failure due to natural wear and tear of any component is not warranted.

OEM factory build-in options (ie: fibre optics, scalers, motors, separators, CVS Kits etc) carry a warranty as recommended by the individual manufacturer.

2 Year Warranted Products

All Belmont intra oral x-rays are warranted for two years from the date of delivery to the original owner. The warranty covers mechanical, electrical and workmanship faults, under normal working conditions. Paintwork and plating are not warranted unless they are reported to Takara Belmont (UK) and proved to be defective at the time of delivery.

Takara Belmont (UK) Ltd are not responsible for any failure caused by misuse of the equipment or damage caused due to poor installation by a non-authorized company.

The warranty policy applies to equipment purchased from and installed by an authorised Takara Belmont (UK) Ltd dealer. Labour will be charged on warranty work after the first year.

Hygiene & Care

In order to maintain & prolong the longevity of your Belmont equipment, we highly recommend using our Hygiene & Care products. (See page 41)

PLEASE NOTE: ISO 13485 Certification

All of our dental products are supplied with a full CE mark and certification to the standard complying with ISO 13485. Any modifications or alterations made after certification may affect your warranty as the responsibility for CE and electrical safety will transfer to the modifier.

| TREATMENT CENTRES (Over-the-Patient delivery system) | With Unit Mounted 921 Light | With Ceiling Mounted D-902C Light | Without Dental Light |
|---|-----------------------------|-----------------------------------|----------------------|
| Eurus S1 'Surgery System' Holder Type (A) | £21,815 | £22,270 | £20,095 |
| Eurus S1 'Surgery System' Holder Type (E) | £24,790 | £25,245 | £23,070 |
| Eurus S1 'Surgery System' CM Rod Type (A) | £23,050 | £23,505 | £21,330 |
| Eurus S1 'Surgery System' CM Rod Type (E) | £26,025 | £26,480 | £24,305 |
| Eurus S1 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System | £22,360 | £22,815 | £20,640 |
| Eurus S1 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System | £25,335 | £25,790 | £23,615 |
| When light sold with package above the saving over recommended retail price is: | Save £1,240 | Save £975 | £0 |

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Chair Mounted Over-Patient Holder & Rod Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation

consists of ~

Chair section features include:

The Eurus S1 has a maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest with more precise positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Integrated service centre
Synchronised axis backrest & seat movement
Integrated foot stalks for chair control
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Autoclavable unit handle cover, handpiece holders
unit tray silicon mat, handpiece rest (Rod Type)
Easy-clean membrane switch panel
Doctor handpiece control touch screen with water low warning indicator
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Electronic foot control for accurate turbine speed control and integrated lifting handle

Unit section (E) Doctor side additional

features include:

Bien-Air MX2 fibre optic micro motor
Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system
2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 180 degree rotation
921 LED dental light with detachable autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Bottled water system with water low warning
LED on cup filler spout casing
Durr spittoon valve



Eurus S1 Flexible System
with Mobile Cart

Flexible Treatment Centres Type (A and E)

with 921 unit-mounted operating light

Optional Right/Left Handed Installation
consists of ~

Chair section features include:

- The Eurus S1 has a maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
- Twin-articulating headrest with more precise positioning and lighter control adjustment
- 2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
- Swing-out armrest to 90 or 135 degrees
- Integrated service centre
- Synchronised axis backrest & seat movement
- Integrated foot stalks for chair control
- Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurus Cabinet Module or Mobile Cart delivery systems A & E (please specify)

'A' version features include:

- Cabinet mounting bracket and arm or height adjustable mobile cart base
- Autoclavable unit handle cover, handpiece holders and unit tray silicon mat
- Easy clean membrane switch panel
- Doctor handpiece control touch screen
- 3-way autoclavable syringe
- 3 Fibre optic & 1 non optic handpiece outlets with individual air & water controls
- First priority handpiece selection
- Electronic foot control for accurate turbine speed control and integrated lifting handle
- Handpiece flush-out
- Handpiece oil mist collector
- Bottled water system

'E' version additional features include:

- Bien-Air MX2 fibre optic micro motor
- Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system
- 2 fibre-optic outlets

Assistant side (A and E)

features include:

- Easily-detached porcelain bowl with 180 degree rotation
- 921 LED dental light with detachable autoclavable handles
- Assistant console with:
 - Extendable console arm safety lockout if arm is lifted
 - Easy-clean membrane chair pre-set controls
 - Detachable HVE & saliva ejector tubings
 - HVE with ball joint for greater ease of operation
 - Automatic bowl-flush & cup filler
 - Bottled water system with water low warning
 - LED on cup filler spout casing
 - Durr spittoon valve

Important:

Maximum loom length between module or cart to the chair is 6.7metres (standard is 2.7metres).

Please read before placing your order:

When ordering a Flexible Treatment Centre, it is imperative that you order the correct item - either a **Mobile Cart or Cabinet Mounted Module.**

If the **incorrect item** has been entered onto your purchase order and delivered to you, the cuspidor and delivery unit will need **returning to Takara Belmont (UK) Ltd** (before unpacking). Two completely new items will need to be manufactured to the correct specification, which will cause delays for your installation at the dental practice, and cost implications may be incurred.

Eurus S1

Eurus S1 Factory Built-in Options



Eurus S1 shown with Luxury Upholstery

for Eurus S1

| | |
|--------|---|
| £0 | NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only) |
| £1,430 | Additional Bien Air MX2 fibre-optic micro motor * |
| £1,705 | Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence |
| £1,598 | Additional NSK NLX-Plus fibre-optic micro motor * |
| £2,145 | Cavitron G-139 scaler |
| £1,465 | NSK Varios 170 non-optic scaler |
| £1,740 | NSK Varios 170 optic scaler |
| £1,530 | Satelec Newtron scaler |
| £1,805 | Satelec Newtron LED scaler |
| £1,020 | Satelec mini LED curing light |
| £2,915 | Durr CAS 1 Combi-Separator with rinse |
| £3,840 | Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount) |
| £5,135 | Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount) |
| £945 | Wireless unit foot controller (when factory fitted) |
| £710 | Monitor mount (if not supplied with camera kit) |
| £410 | Additional 3-way autoclavable syringe (assistant side) |
| £198 | Additional easily detachable HV ejector hose (including ball joint for greater ease of operation) |
| £185 | Large tray (double size for Rod Type) |
| £90 | Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm) |
| £235 | Flushing Tray A (Holder Type) |
| £235 | Flushing Tray B (Rod Type) |
| £255 | Annual service Kit 7 |

* Additional micro motors must be the same model as the primary micro motor.



New Flushing Tray - suitable for all treatment centres

| TREATMENT CENTRES (Below-the-Patient delivery system) | With Unit Mounted 921 Light | With Ceiling Mounted D-902C Light | Without Dental Light |
|---|-----------------------------|-----------------------------------|----------------------|
| Eurus S6 'Surgery System' Holder Type (A) | £27,620 | £28,075 | £25,900 |
| Eurus S6 'Surgery System' Holder Type (E) | £30,595 | £31,050 | £28,875 |
| Eurus S6 'Vac-Pac System' Holder Type (A) | £22,915 | £23,370 | £21,195 |
| Eurus S6 'Vac-Pac System' Holder Type (E) | £25,890 | £26,345 | £24,170 |
| Eurus S6 'Flexible System' (A) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System | £26,920 | £27,375 | £25,200 |
| Eurus S6 'Flexible System' (E) Dental Chair and Cuspidor with Eurus Cabinet Mounted Module or Mobile Cart Delivery System | £29,895 | £30,350 | £28,175 |
| When light sold with package above the saving over recommended retail price is: | Save £1,240 | Save £975 | £0 |

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Surgery System (A and E) with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Taller backrest cushion with additional lumbar support for greater patient comfort

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system

Height adjustable unit tray with pneumatic release and lock button

Powerful intuitive touch-screen with advanced instrument control

Autoclavable unit handle cover, handpiece holders unit tray silicon mat

Easy-clean membrane switch panel

3-way autoclavable syringe

3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control

First priority handpiece selection

Handpiece water line flush-out system

Handpiece oil mist collector

Wireless electronic foot control for accurate instrument operation with integrated lifting handle

Bottle water system

Unit section (E) Doctor side additional features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with auto-reverse, Endo control and LCD display, with optional up-grade to WaveOne endodontic system 2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with 90 degree rotation

Cuspidor with 90 degree rotation away from the chair for accessing the patient from the left side 921 LED dental light with detachable autoclavable handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the cup filler spout and touch-screen

Durr spittoon valve

'Vac-Pac System' (A and E)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

Assistant's console

Easy-access detachable filter

Autoclavable HVE & saliva ejector

Funnel connection hose



Eurus S6 enables greater workflow freedom

Eurus S6 Treatment Centres with folding and extending legrest



Flexible Treatment Centres Type (A and E)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 165kg
(IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise
positioning and lighter control adjustment

2 pre-set positions, last position memory (LP) &
auto-return plus 6 operator configurations

Chair armrest folds away 180 degrees towards the
back of the chair for easy infirm patient entry

Synchronised axis backrest & seat movement

Taller backrest cushion with additional lumbar
support for greater patient comfort

Integrated foot stalks for chair control

Integrated service centre

Seamless upholstery as standard

'Flexible System A & E'

Flexible System with Eurus Cabinet Module or
Mobile Cart delivery systems A & E (please specify)

'A' version features include:

Cabinet mounting bracket and arm or height
adjustable mobile cart base

Autoclavable unit handle cover, handpiece holders
and unit tray silicon mat

Easy clean membrane switch panel

Doctor handpiece control touch screen

3-way autoclavable syringe

3 Fibre optic & 1 non optic handpiece outlets with
individual air & water controls

First priority handpiece selection

Wireless electronic foot control for accurate
instrument operation with integrated lifting handle

Handpiece flush-out

Handpiece oil mist collector

Bottled water system

'E' version additional features include:

Bien-Air MX2 fibre optic micro motor

Pre-set micro motor & torque control with
auto-reverse, Endo control and LCD display, with
optional up-grade to WaveOne endodontic system

2 fibre-optic outlets

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with
90 degree rotation

Cuspidor with 90 degree rotation away from the
chair for accessing the patient from the left side

921 LED dental light with detachable autoclavable
handles

Assistant console with:

Extendable console arm safety lockout if arm is lifted

Easy-clean membrane chair pre-set controls

Detachable HVE & saliva ejector tubings

2 x HVE with ball joint for greater ease of operation

Automatic bowl-flush & cup filler

Chair movement lock indicators located on the
cup filler spout and touch-screen

Durr spittoon valve

Important:

**Maximum loom length between module
or cart to the chair is 6.7metres
(standard is 2.7metres).**

Please read before placing your order:

When ordering a Flexible Treatment Centre, it is
imperative that you order the correct item - either
a **Mobile Cart or Cabinet Mounted Module**.

If the **incorrect item** has been entered onto your
purchase order and delivered to you, the cuspidor
and delivery unit will need **returning to Takara**

Belmont (UK) Ltd (before unpacking). Two
completely new items will need to be manufactured
to the correct specification, which will cause delays
for your installation at the dental practice, and cost
implications may be incurred.

Eurus S6



Eurus S6 shown with Luxury Upholstery

Eurus S6 Factory Built-in Options



for Eurus S6

| | |
|--------|--|
| £0 | NSK NLX-Plus micro motor (in place of Bien Air MX2 micro motor on 'E' version only) |
| £1,430 | Additional Bien Air MX2 fibre-optic micro motor * |
| £1,705 | Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence |
| £1,598 | Additional NSK NLX-Plus fibre-optic micro motor * |
| £2,145 | Cavitron G-139 scaler |
| £1,465 | NSK Varios 170 non-optic scaler |
| £1,740 | NSK Varios 170 optic scaler |
| £1,530 | Satelec Newtron scaler |
| £1,805 | Satelec Newtron LED scaler |
| £1,020 | Satelec mini LED curing light |
| £2,915 | Durr CAS 1 Combi-Separator with rinse |
| £3,840 | Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount) |
| £5,135 | Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount) |
| £710 | Monitor mount (if not supplied with camera kit) |
| £410 | Additional 3-way autoclavable syringe (assistant side) |
| £90 | Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm) |
| £235 | Flushing Tray A (Holder Type) |
| £255 | Annual Service Kit 7 |

* Additional micro motors must be the same model as the primary micro motor.

| TREATMENT CENTRES (Below-the-Patient delivery system) | With Unit Mounted 921LR Light | With Ceiling Mounted D-902C Light | Without Dental Light |
|--|----------------------------------|--------------------------------------|----------------------------|
| Eurus S8 'Surgery System' Holder Type (A) | £23,845 | £24,300 | £22,125 |
| Eurus S8 'Surgery System' Holder Type (E) | £26,815 | £27,270 | £25,095 |
| Eurus S8 'Vac-Pac System' Holder Type (A) | £21,565 | £22,020 | £19,845 |
| Eurus S8 'Vac-Pac System' Holder Type (E) | £24,540 | £24,995 | £22,820 |
| Eurus S8 'Surgery System' CM Rod Swing (A) | £26,980 | £27,435 | £25,260 |
| Eurus S8 'Surgery System' CM Rod Swing (E) | £29,950 | £30,405 | £28,230 |
| Eurus S8 'Vac-Pac System' CM Rod Swing (A) | £23,995 | £24,450 | £22,275 |
| Eurus S8 'Vac-Pac System' CM Rod Swing (E) | £26,970 | £27,425 | £25,250 |
| When light sold with package above the saving over recommended retail price is: | Save £1,240 | Save £975 | £0 |

Please state light model number when ordering treatment centre

Please see swatch brochure for upholstery options

Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg
(IEC60601-1 / ISO 6875)

Twin-articulating headrest with more precise
positioning and lighter control adjustment
2 pre-set positions, last position memory (LP) &
auto-return plus 6 operator configurations
Swing-out armrest to 90 or 135 degrees
Synchronised axis backrest & seat movement
Integrated foot stalk or combined with unit foot
controller (Holder Type)
Integrated foot stalks for chair control (Rod Type)
Integrated service centre
Seamless upholstery as standard

Unit section (A) Doctor side: features include:

Below-the-patient swing-arm delivery system
Height adjustable unit tray with 3 settings (Holder
Type)
Powerful intuitive touch-screen with advanced
instrument control
Autoclavable unit handle cover, handpiece holders
& tray silicon mat / Rod Type handpiece place mat
Easy-clean membrane switch panel
3-way autoclavable syringe
3 Fibre-optic & 1 non-optic handpiece outlets with
individual air & water control
First priority handpiece selection
Handpiece water line flush-out system
Handpiece oil mist collector
Wireless electronic foot control for accurate
instrument operation with integrated lifting handle
plus chair manual control movements (Holder Type)
Bottle water system

Unit section (E) Doctor side additional

features include:

NSK NLX-Plus fibre optic micro motor
Pre-set micro motor & torque control with
auto-reverse & touch-screen display
2 fibre-optic outlets

WaveOne is available as an additional option
when the Bien Air MX2 micro motor is fitted

Assistant side (A and E)

features include:

Easily-detached porcelain bowl with
90 degree rotation
Cuspidor with 90 degree rotation away from the
chair for accessing the patient from the left and
right sides
921LR LED dental light with detachable
autoclavable handles
Assistant console with:
Extendable console arm safety lockout if arm is lifted
Easy-clean membrane chair pre-set controls
Detachable HVE & saliva ejector tubings
2 x HVE with ball joint for greater ease of operation
Automatic bowl-flush & cup filler
Chair movement lock indicators located on the
cup filler spout and touch-screen
Durr spittoon valve



Base mounted ambidextrous
Holder-Type (A and E) treatment centre

Eurus S8 shown with Luxury Upholstery

Surgery System (A and E)

(Ambidextrous)

with 921 unit-mounted operating light

consists of ~

Chair section features include:

Maximum safe working patient weight of 200kg
(IEC60601-1 / ISO 6875)

- Twin-articulating headrest with more precise positioning and lighter control adjustment
- 2 pre-set positions, last position memory (LP) & auto-return plus 6 operator configurations
- Swing-out armrest to 90 or 135 degrees
- Synchronised axis backrest & seat movement
- Integrated foot stalk or combined with unit foot controller (Holder Type)
- Integrated foot stalks for chair control (Rod Type)
- Integrated service centre
- Seamless upholstery as standard

Unit section (A) Doctor side: features include:

- Below-the-patient swing-arm delivery system
- Height adjustable unit tray with 3 settings (Holder Type)
- Powerful intuitive touch-screen with advanced instrument control
- Autoclavable unit handle cover, handpiece holders & tray silicon mat / Rod Type handpiece place mat
- Easy-clean membrane switch panel
- 3-way autoclavable syringe
- 3 Fibre-optic & 1 non-optic handpiece outlets with individual air & water control
- First priority handpiece selection
- Handpiece water line flush-out system
- Handpiece oil mist collector
- Wireless electronic foot control for accurate instrument operation with integrated lifting handle plus chair manual control movements (Holder Type)
- Bottle water system

Unit section (E) Doctor side additional

features include:

- NSK NLX-Plus fibre optic micro motor
- Pre-set micro motor & torque control with auto-reverse & touch-screen display
- 2 fibre-optic outlets

WaveOne is available as an additional option when the Bien Air MX2 micro motor is fitted

'Vac-Pac System' (A and E)

(Ambidextrous)

includes the same as 'Surgery System' with a Vac-Pac system (in lieu of a cuspidor) features include:

- Assistant console with:
 - Extendable console arm safety lockout if arm is lifted
 - Easy-access detachable filter
 - Autoclavable HVE & saliva ejector
 - Funnel connection hose
 - Chair control auto & manual touch pad
- Assistant tray

Eurus S8

Chair mounted under patient ambidextrous Rod Swing (A and E) Treatment Centre shown with optional large tray



Eurus S8 shown with Luxury Upholstery

Eurus S8 Factory Built-in Options



for Eurus S8

| | |
|--------|--|
| £0 | Bien Air MX2 micro motor (in place of NSK NLX-Plus micro motor on 'E' version only) |
| £1,430 | Additional Bien Air MX2 fibre-optic micro motor * |
| £1,705 | Upgrade to Bien Air WaveOne Contra Angle, plus one-off activation licence |
| £1,598 | Additional NSK NLX-Plus fibre-optic micro motor * |
| £2,145 | Cavitron G-139 scaler |
| £1,465 | NSK Varios 170 non-optic scaler |
| £1,740 | NSK Varios 170 optic scaler |
| £1,530 | Satelec Newtron scaler |
| £1,805 | Satelec Newtron LED scaler |
| £1,020 | Satelec mini LED curing light |
| £2,915 | Durr CAS 1 Combi-Separator with rinse |
| £3,840 | Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount) |
| £5,135 | Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount) |
| £710 | Monitor mount (if not supplied with camera kit) |
| £410 | Additional 3-way autoclavable syringe (assistant side) |
| £185 | Large tray (double size for Rod Type) |
| £90 | Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm) |
| £235 | Flushing Tray A (Holder Type) |
| £235 | Flushing Tray B (Rod Type) |
| £255 | Annual Service Kit 7 |

* Additional micro motors must be the same model as the primary micro motor.

| TREATMENT CENTRES Below-the-Patient delivery system and flexible options | With Unit Mounted 320LR Light | With Ceiling Mounted 302C Light | Without Dental Light |
|--|-------------------------------|---------------------------------|----------------------|
| Voyager III 'Surgery System' | £15,360 | £15,760 | £14,020 |
| Voyager III 'SP' Chair with ambidextrous cuspidor with Voyager cabinet mounted module or mobile cart delivery system | £15,360 | £15,760 | £14,020 |
| When light sold with package above the saving over recommended retail price is: | Save £855 | Save £870 | £0 |

| TREATMENT CENTRES Below-the-Patient delivery system and flexible options | With Unit Mounted 920LR Light | With Ceiling Mounted D-902C Light | Without Dental Light |
|---|-------------------------------|-----------------------------------|----------------------|
| Voyager III 'Surgery System' | £15,935 | £16,745 | £14,020 |
| Voyager III 'SP' As above | £15,935 | £16,745 | £14,020 |
| When light sold with package above the saving over recommended retail price is: | Save £860 | Save £425 | £0 |

| | |
|---|---------|
| Voyager III 'SP' Chair with ambidextrous cuspidor (without operating light) | £10,750 |
|---|---------|

Please state light model number when ordering treatment centre

MODULE & MOBILE CART OPTIONS

| | |
|---|--------|
| Voyager 'M' & 'MC' Module and Mobile Cart delivery system | £3,270 |
|---|--------|

'Surgery System' (Ambidextrous)

with 320LR unit-mounted Operating Light

consists of ~

Chair section features include:

The Voyager III has a maximum safe working patient weight of 135kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrests (left & right)
Integral service centre with air & water filter regulators
8-way foot control

Base-mounted swing-around delivery system features include:

3 handpiece outlets with individual air & water control
3-way autoclavable syringe
Table safety lock device
Quick flush-out system
Bottled water system
Stainless steel instrument tray
Non-slip mat

Chair-mounted rotatable cuspidor unit

features include:

Porcelain cuspidor bowl
Manual cup-filler & bowl-flush
1 detachable high volume evacuator
1 detachable saliva ejector
Easily-accessed & detached solids collector
Durr spittoon valve

'SP' (Ambidextrous)

Voyager III Dental Chair with cuspidor

consists of ~

Chair section features include:

The Voyager III has a maximum safe working patient weight of 135kg (IEC60601-1 / ISO 6875)
Twin-articulating headrest
2 preset positions, Last Position Memory (LP) & Auto-Return
Swing-out armrests (left & right)
Integral service centre with air & water filter regulators
8-way foot control

Chair-mounted rotatable cuspidor unit features include:


Porcelain cuspidor bowl
Manual cup-filler & bowl-flush
1 detachable high volume evacuator
1 detachable saliva ejector
Easily-accessed & detached solids collector
Durr spittoon valve

Voyager III Factory Built-in Options

for Voyager

| | |
|--------|--|
| £2,505 | Bien Air Optima MCX LED (inc. standard fibre-optic micro motor, turbine hose & touch pad control) |
| £2,505 | NSK NLX-Nano fibre-optic pack (inc. NLX-Nano fibre-optic micro motor, turbine hose & touch pad control) |
| £3,055 | NSK NLX-Plus multi-pad endo pack (inc. NLX-Plus fibre-optic micro motor, turbine hose & multi-pad Endo control) |
| £870 | Bien Air fibre-optic kit (inc. 1 turbine hose) |
| £298 | Additional fibre-optic outlet |
| £2,245 | Cavitron G-139 scaler |
| £1,465 | NSK Varios 170 non-optic basic scaler |
| £1,740 | NSK Varios 170 optic scaler (with touch pad control) |
| £1,565 | Satelec Newtron scaler |
| £1,850 | Satelec Newtron LED scaler |
| £1,215 | Satelec mini LED curing light |
| £440 | Additional 3-way autoclavable syringe on Assistant's side inc. bottled water system |
| £0 | Standard (10x10) Module Arm (alternative to Voyager Arm System) |
| £515 | Touch pad assembly for chair control fitted to Operator's table |
| £830 | Monitor mount (if not supplied with camera kit) |
| £3,840 | Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount) |
| £5,135 | Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount) |
| £85 | Extra-long umbilical (up to 4 metres for Mobile Cart only) |
| £90 | Tissue box holder (suits 72 wipe box H108mm, W138mm, D31mm) |
| £120 | Floor box fitted to cart base |
| £235 | Flushing Tray A (Holder Type) |
| £245 | Annual Service Kit 1 (for Voyager I, II & III Treatment Centres) |
| £158 | Annual Service Kit 2 (for Voyager Module and Cart) |

Did you know you can have a private viewing?



Our showrooms are open during office hours, weekends & evenings, whilst our 3D showroom experience is available from the comfort of your home.

You can access our showroom 24/7 from wherever you are, simply go to belmontdental.co.uk/showroom

Flexible Delivery Options

MODULE & MOBILE CART OPTIONS

| | |
|---|---------------|
| Voyager 'M' & 'MC' Module and Mobile Cart delivery system | £3,270 |
| Clesta II 'M' & 'MC' (A) Module and Mobile Cart delivery system | £5,195 |
| Clesta II 'M' & 'MC' (E) Module and Mobile Cart delivery system | £8,215 |

Universally compatible
Flexible Options



Voyager 'M'

Module delivery system features include:

- Mounting arms & bracket
- 3 handpiece outlets with individual air & water control
- 3-way autoclavable syringe
- Bottled water system
- Stainless steel instrument tray
- Non-slip mat
- Service centre with air filter regulator

Voyager 'MC'

Mobile Cart delivery system includes the same as Voyager 'M' but with height-adjustable mobile base

Clesta II 'M' (A)

Module delivery system features include:

- Mounting arms & bracket
- 3-way autoclavable syringe
- 3 handpiece outlets with individual air & water control
- Handpiece flush-out system
- Bottled water system
- Service centre with air filter regulator

Clesta II 'MC' (A)

Mobile Cart delivery system includes the same as Clesta II 'M' (A) but with height-adjustable mobile base

Clesta II 'M' (E)

Module delivery system includes the same as the 'A' type but with the following additional features:

- NSK NLX-Plus fibre-optic micro motor
- Preset micro motor & digital display
- 2 fibre-optic outlets
- First Priority handpiece selection
- Handpiece oil mist collector
- Additional chair preset positioning

Clesta II 'MC' (E)

Mobile Cart delivery system includes the same as Clesta II 'M' (E) but with height-adjustable mobile base

Clesta II

'Flexible System'

Module Delivery System



Mobile Cart Delivery System

Clesta II Module & Cart **Factory Built-in Options**



for Clesta II

| | |
|---------------|--|
| £0 | Bien Air MX2 micro motor (in place of NSK NLX-Plus micro motor on 'E' version only) |
| £3,615 | Bien Air MX2 Optima INT Endo pack (inc. standard fibre-optic micro motor, turbine hose & Optima INT Endo control) |
| £2,480 | NSK NLX-Nano fibre-optic pack (inc. NLX-Nano fibre-optic micro motor & turbine hose) |
| £3,055 | NSK NLX-Plus multi-pad Endo pack (inc. NLX-Plus fibre-optic micro motor, turbine hose & multi-pad Endo control) |
| £2,130 | Additional Bien Air MX2 fibre-optic micro motor * |
| £1,120 | Additional NSK NLX-Nano fibre-optic micro motor * |
| £1,598 | Additional NSK NLX-Plus fibre-optic micro motor * |
| £810 | Bien Air fibre-optic kit (inc. 1 turbine hose) |
| £298 | Additional fibre-optic outlet |
| £2,145 | Cavitron G-139 scaler |
| £1,465 | NSK Varios 170 non-optic scaler |
| £1,740 | NSK Varios 170 optic scaler |
| £1,530 | Satelec Newtron scaler |
| £1,805 | Satelec Newtron LED scaler |
| £1,020 | Satelec mini LED curing light |
| £3,840 | Satelec 617 camera (inc. USB2 connection, freeze frame & monitor mount) |
| £5,135 | Satelec Soprocure camera (inc. USB2 connection, freeze frame & monitor mount) |
| £170 | Handpiece oil mist collector (Flexible 'A' version only) |
| £85 | Extra-long umbilical (up to 4 metres for Mobile Cart only) |
| £120 | Floor box fitted to cart base |
| £175 | Annual Service Kit 3 (for Clesta I & II Module and Cart) |

Dental Chairs

WITH FOLDING AND EXTENDING LEGREST

| | |
|---|----------------|
| Eurus Swivel Chair Standalone chair <i>with</i> Twin-articulating headrest & cable free base rotation 105 degrees (L & R) available in all seamless and luxury upholstery | £9,835 |
| Eurus Swivel PH Chair <i>with</i> Power headrest (special order) available in all seamless and luxury upholstery | £10,905 |
| 050 Clair Chair <i>with</i> Twin-articulating headrest available in all seamless and luxury upholstery | £9,355 |
| 055 Clair Chair <i>with</i> Power headrest available in all seamless and luxury upholstery | £10,345 |

WITH FIXED LEGREST

| | |
|--|---------------|
| Clesta II Chair Standalone chair | £7,035 |
|--|---------------|

Folding & Fixed Legrest Dental Chairs



Eurus Swivel Chair

features include:

Standalone chair
 Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)
 Cable free base rotation 105 degrees (L & R)
 Twin-articulating headrest with more precise positioning and lighter adjustment
 Synchronised axis backrest, seat & leg-rest movement
 Integrated foot stalks for chair control with 2 pre-set positions, Last Position Memory (LP) & Auto-return
 Automatic extendable leg-rest with manual override & auto safety stop
 Foldaway armrests
 Seamless upholstery as standard

Eurus Swivel PH Chair (special order)

features include the same as above but with a Power headrest

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

050 Clair Chair

features include:

Maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
 Twin-articulating headrest
 2 preset positions, Last Position Memory (LP) & Auto-Return
 Automatic extendable legrest with manual override & auto-safety stop
 Foldaway, detachable armrests
 Seamless upholstery

055 Clair Chair

features include the same as 050 Clair Chair but with a Power headrest

Maximum safe working patient weight of 165kg (IEC60601-1 / ISO 6875)

Clesta II Chair

features include:

Maximum safe working patient weight of 200kg (IEC60601-1 / ISO 6875)
 Twin-articulating headrest
 2 preset positions, Last Position Memory (LP) & Auto-Return
 Swing-out armrest
 Synchronised axis backrest
 8-way foot control
 Seamless upholstery
 Standalone chair

All chairs above in this column can have a light post and bracket fitted (see page 35 for details)

Intra-Oral X-Ray Units

097 BELRAY II

| | |
|---|---------------|
| 097WK-300 Belray II wall-mounted, with 300mm primary arm, wall plate, scissor-type arm | £3,465 |
| 097WK-800 Belray II wall-mounted, with 800mm primary arm, wall plate, scissor-type arm | £3,465 |
| 097WK-1000 Belray II wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm | £3,465 |

Belmont Touch PHOT-X IIS DC LCD X-RAY

| | |
|---|---------------|
| 505WK-300 wall-mounted, with 300mm primary arm, wall plate, scissor-type arm | £4,125 |
| 505WK-800 wall-mounted, with 800mm primary arm, wall plate, scissor-type arm | £4,125 |
| 505WK-1000 wall-mounted, with 1,000mm primary arm, wall plate, scissor-type arm | £4,125 |
| 505FM mobile base | £4,535 |
| 505CK ceiling-mounted, with ceiling-mount pole & bracket | £4,975 |
| 505FK floor-mounted base | £3,855 |

097 Belray II & Phot-X IIS DC LCD Intra-Oral X-Ray Units & optional extras



097 Belray II

features include:

- 70 KVP
- 0.7mm focal spot
- Digital time display
- Digital imaging compatible (min. exposure 0.02 sec)
- Manual override for total flexibility
- Touch control for 2 film settings & digital imaging
- Lockable touch pad control
- Optional cone length & tube current
- Rectangular Collimator
- Hand exposure switch assembly

models include:

097WK-300

097WK-800

097WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach) of 300mm (1,350*) / 800mm (1,850*) / 1,000mm (2,050*)

Optional extras for 097 Belray II

£325 Long cone for paralleling technique

Optional extras for Phot-X IIS

£325 Long cone for paralleling technique

£345 Operating Light bracket for 505CK

Phot-X IIS DC LCD X-Ray

features include:

- Touchscreen control
- Easy-clean glass screen
- Safety-conscious technologies
- USB port
- Instant Error Code reports
- Lockable keypad
- Adult/adolescent/child presets
- 60 or 70 KVP options
- 0.4mm focal spot
- Digital time display
- Radiation dose display
- Compatible for digital imaging (min. exposure 0.01 sec)
- Manual override for total flexibility
- Optional cone length, tube current & tube voltage
- Rectangular Collimator
- Hand exposure switch assembly

models include:

505WK-300

505WK-800

505WK-1000

Wall-mounted, wall plate, scissor-type arm

*** OPTIONS:** primary arm length (*Max. reach) of 300mm (1,420*) / 800mm (1,920*) / 1,000mm (2,120*)

505FM Mobile base

505FK Floor-mounted base

505CK Ceiling-mounted, with ceiling-mount pole & bracket

*** OPTIONS:**

| Ceiling Height (mm) | Pole Length (mm) |
|---------------------|------------------|
| 2,700 | 1,070 |
| 3,000 | 1,370 |
| 3,200 | 1,570 |

***Please specify when ordering**

Operating Lights

300 LED SERIES MODELS

| | |
|---|---------------|
| 301U LED post-mounted | £1,805 |
| 302C LED ceiling-mounted | £2,610 |
| 320LR LED unit-mounted | £2,195 |
| 320LR LED + Kit unit-mounted light plus fitting kit 1P051PA0-KIT (to replace the 048 & 720S LR lights for Voyager II/III only) | £2,395 |

900 LED SERIES MODELS

| | |
|--|---------------|
| 901 LED post-mounted | £2,835 |
| D-902C LED ceiling-mounted, suitable for a maximum ceiling height 305cm | £3,150 |
| D-902CX LED ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray | £2,895 |
| 905 LED ceiling-track - please specify pole length required | £4,100 |
| 920LR LED chair-mounted, only compatible with Voyager III & III 'SP' | £2,775 |
| 921 LED unit-mounted, only compatible with Eurus S1 & S6 | £2,960 |
| 921LR LED unit-mounted, only compatible with Eurus S8 | £2,960 |

LED Operating Lights & optional extras



300 LED Series

features include:

- White natural beam, 5000 Kelvin (5 LEDs)
- Touchless sensor switch
- Composite Curing Safe Mode
- Triple-axis head movement
- Auto ON/OFF with preset chair positions*

models include:

- 301U LED** Post-mounted
- 302C LED** Ceiling-mounted
- 320LR LED** Unit-mounted, *only compatible with Voyager III 'Surgery System' & Voyager III 'SP'*

* Applies to unit-mounted versions only.

900 LED Series

features include:

- White natural beam, 5000 Kelvin (10 LEDs)
- Touchless sensor switch
- Composite Curing Safe Mode
- Triple-axis head movement
- Auto ON/OFF with preset chair positions

models include:

- 901 LED** Post-mounted
- D-902C LED** Ceiling-mounted, suitable for a maximum ceiling height 305cm
- D-902CX LED** Ceiling-mounted, to be used with Belmont ceiling-mounted X-Ray unit
- 905 LED** Ceiling-track with pole:

* OPTIONS:

| Ceiling Height (mm) | Pole Length (mm) |
|---------------------|------------------|
| 2,380 | 380 |
| 2,480 | 480 |
| 2,530 | 530 |
| 2,680 | 680 |
| 2,840 | 840 |
| 3,000 | 1,000 |

920LR LED Chair-mounted, only compatible with Voyager III and III 'SP'

921 LED Unit-mounted, only compatible with Eurus S1 & S6

921LR LED Unit-mounted, only compatible with Eurus S8

***Please specify when ordering**

Optional extras for Belmont Operating Lights

| | |
|-------------|--|
| £560 | Chair-mounted light-post & bracket for 050 Clair and Clesta II standalone chairs (Part No. UK-149) |
| £205 | 1,500 mm ceiling-pole in lieu of 1,000 mm pole (suitable for max. ceiling height 3,550 mm) (Part No. UK-068) |
| £130 | Patient Mirror for 900 LED Series (Part No. 1G02WMB0) |

Dental Stools

SADDLE STOOL

| | |
|---------------------------------|-------------|
| Saddle Stool Gas Lift | £620 |
|---------------------------------|-------------|

FOR OPERATORS

| | |
|-------------------------|-------------|
| ic-D Gas Lift | £435 |
|-------------------------|-------------|

| | |
|--------------------------|-------------|
| ic-DF Gas Lift | £460 |
|--------------------------|-------------|

| | |
|----------------------------|-------------|
| ic-Tilt Gas Lift | £605 |
|----------------------------|-------------|

FOR ASSISTANTS

| | |
|---------------------------|-------------|
| ic-DSA Gas Lift | £470 |
|---------------------------|-------------|

| | |
|----------------------------|-------------|
| ic-DSAF Gas Lift | £565 |
|----------------------------|-------------|

| | |
|----------------------------------|-------------|
| ic-DSSA Swing-Arm type | £650 |
|----------------------------------|-------------|

Operators' & Assistants' Dental Stools & optional extras

Saddle Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Adjustable seat tilt angle

Seat height 580 - 780mm max.

ic-D Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height

Seat height 490 - 620mm max.

ic-DF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle

Seat height 500 - 630mm max.

ic-Tilt Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable seat & backrest angle

Seat height 500 - 630mm max.

ic-DSA Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Footrest ring

Seat height 530 - 730mm max.

ic-DSAF Stool

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Free floating and lockable backrest angle
Footrest ring

Seat height 540 - 740mm max.

ic-DSSA Stool

Swing-Arm type

Five-point polished aluminium base
Smooth running castors
Gas lift adjustable seat height
Independent backrest height adjustment
Backrest horizontal forward and back movement
Footrest ring

Seat height 590 - 730mm max.

* All stools benefit from:

Easy clean upholstery
Double-density cushioning

Optional extras for Operating Stools

| | |
|------------|-------------------------------------|
| £55 | Footring for ic-Tilt , ic-DF & ic-D |
|------------|-------------------------------------|

Please see swatch brochure for upholstery options

Upholstery

Fabric range options for Treatment Centres, Dental Chairs and Dental Stools

SEAMLESS UPHOLSTERY

| | |
|-------------------------------------|-------------|
| Standard colours | £0 |
| Ultrasoft colours | £420 |
| Ultrasoft Pro colours | £420 |
| Ultrasoft Pearlished colours | £420 |

LUXURY UPHOLSTERY

| | |
|---|-------------|
| Ultrasoft colours Not available on Voyager III treatment centres | £710 |
| Ultrasoft Pro colours Not available on Voyager III treatment centres | £710 |
| Ultrasoft Pearlished colours Not available on Voyager III treatment centres | £710 |

DENTAL STOOLS

| | |
|-------------------------------------|------------|
| Standard colours | £0 |
| Ultrasoft colours | £95 |
| Ultrasoft Pro colours | £95 |
| Ultrasoft Pearlished colours | £95 |

IMPORTANT NOTE:

Cancelled or amended orders of Ultrasoft, Ultrasoft Pro or Ultrasoft Pearlished upholstery will incur the following fee: £500 per Luxury chair, £250 per Seamless chair, and £50 per stool.

Seamless Upholstery

Available in a choice of four fabrics



Eurus S1 Chair shown with Seamless Upholstery

Luxury Upholstery

Available in a choice of three fabrics



Eurus Swivel Chair shown with Luxury Upholstery

Hygiene & Care

For Cleaning Standards Beyond the Standard

For cleaning and maintenance of your treatment centre we highly recommend the following products:

TREATMENT CENTRE MAINTENANCE

| | |
|---|------------|
| B100 (500ml) Upholstery cleaning and care | £38 |
| B300 (1 litre) Surface disinfectant | £32 |
| B500 DUWL Starter Kit - 2 periodic shock tests, 2 periodic shock treatments and 20 daily maintenance sachets | £54 |

DENTAL UNIT WATERLINE PROTECTION

| | |
|---|------------|
| B700 (100 sachets) DUWL daily maintenance | £98 |
| B900 (5 x 200ml) DUWL periodic shock and test treatment | £78 |

Effective against all viruses, including COVID-19

Designed for easy use

Guarantees safety during application

Approved by the Commission for Hospital Hygiene & Infectious Disease Prevention at the Robert Koch Institute (RKI)

We're offering a **free Intro Kit** with every new treatment centre ordered. This includes a **B100**, **B300** and **B500** as an introduction to our hygiene & care products.



Dental Equipment Service Kits

ANNUAL SERVICE KIT 1

for : Voyager I, II & III Treatment Centres

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK001 | SK for Auto-select block 3-handpieces (1P055WA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK004 x2 | SK for Large/main master valve (1P055ZA0) |
| Y3072 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| DAQR51A0 | Oil Separator Filter for Voyager/Voyager-II L Unit |
| FJOA09A0 | Element for Air Filter Regulator |
| Y8136 | CWS Manifold Washer |
| FJOA08A0 | Element for Water Filter Regulator |
| 1A025DB0 | Drain Strainer |
| DU-7560-500-22 | Diaphragm for Isolation Valve (7560-500-60) |
| £245 | |

ANNUAL SERVICE KIT 2

for : Voyager Module and Cart

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK001 | SK for Auto-select block 3-handpieces (1P055WA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK004 | SK for Large/main master valve (1P055ZA0) |
| Y3072 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| DAQR51A0 | Oil Separator Filter for Voyager/Voyager-II L Unit |
| Y8136 | CWS Manifold Washer |
| 5726-04 | Filter for Voyager-M Air Regulator (B07-201-A1KG) |
| £158 | |

ANNUAL SERVICE KIT 3

for : Clesta I & II Module and Cart

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK002 | SK for Auto-select block 4-handpieces (1P055XA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK006 | SK for Small master valve (1P0561B0) |
| Y3072 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| FJOA09A0 | Element for air filter regulator |
| Y8136 | CWS Manifold Washer |
| FARA71B0 | Sponge for Oil Mist Collector |
| £175 | |

ANNUAL SERVICE KIT 4

for : Clesta II and Compass Treatment Centres

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK002 | SK for Auto-select block 4- handpieces (1P055XA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK004 x2 | SK for Large/main master valve (1P055ZA0) |
| Y3072 x2 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| FJOA09A0 | Element for Air Filter Regulator |
| FJOA08A0 | Water Filter Element |
| Y8136 | CWS Manifold Washer |
| DU-7560-500-22 | Diaphragm for Isolation Valve (7560-500-60) |
| 1A025DB0 | Drain Strainer |
| FARA71B0 | Sponge for Oil Mist Collector |
| £280 | |

Dental Equipment Service Kits

ANNUAL SERVICE KIT 5

for : SP-Cleo and Cleo II Treatment Centres

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK002 | SK for Auto-select block 4- handpieces (1P055XA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK004 x2 | SK for Large/main master valve (1P055ZA0) |
| Y3072 x2 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| FJOA09A0 | Element for Air Filter Regulator |
| FJOA08A0 | Water Filter Element |
| Y8136 | CWS Manifold Washer |
| DU-7560-500-22 | Diaphragm for Isolation Valve (7560-500-60) |
| 1A025DB0 | Drain Strainer |
| FARA71B0 | Sponge for Oil Mist Collector |
| SK006 x2 | SK for small master valve (1P0561B0) |
| £325 | |

ANNUAL SERVICE KIT 6

for : Clesta I Treatment Centres

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| SK002 | SK for Auto-select block 4- handpieces (1P055XA0) |
| SK003-A | SK for Water relay valve (1P06BWA0) |
| SK004 x2 | SK for Large/main master valve (1P055ZA0) |
| Y3072 x2 | Deluxe Autoclavable Valve Core Style Syringe - Syringe Valve O-Ring Kit |
| Y2230 | O-Ring, Buna-n - .042 x 0050 (Pack of 12) |
| FJOA09A0 | Element for Air Filter Regulator |
| FJOA08A0 | Water Filter Element |
| SK005 | SK for Bowl Flush & Cup filler Block (1P0560B0) |
| Y8136 | CWS Manifold Washer |
| DU-7560-500-22 | Diaphragm for Isolation Valve (7560-500-60) |
| 1A025DB0 | Drain Strainer |
| £325 | |

ANNUAL SERVICE KIT 7

for : Eurus S1, S6 & S8 Treatment Centres

| <i>consists of ~</i> | |
|----------------------|---|
| Part Number | Description |
| 1P055ZA0 x2 | Service Kit Water/Air Main Master Valve |
| 1E09J6A0 | PN8136 Washer for Water Bottle |
| 1E097DA0 | Water Filter Element |
| 1E07K9A0 | Air Filter Element 40 |
| 1A0GCEA0 x2 | Solid Collector Filter - Cone Shape |
| 1P07BJA0 | EURUS Water Relay Valve Kit |
| 1A0PJKA0 | Diaphragm for Air Valve |
| 1A0PURA0 | Sponge for Oil Mist |
| 1A025DB0 | Drain Strainer |
| 1E09J4A0 | PN3072 Valve Assy for DCI Syringe |
| 1E09J5A0 | PN2230 O-Ring Kit |
| £255 | |



Scan to download PDF price list:



SALES

LONDON & SOUTH EAST

Eugene O'Malley **T:** 0771 097 9668

MIDLANDS, SOUTH WEST, NORTHERN & SOUTHERN IRELAND

Barry Hunt **T:** 0771 526 3587

NORTHERN ENGLAND & SCOTLAND

Mark Harris **T:** 0786 053 5968

London Head Office & Showroom

Belmont House, One St. Andrews Way, Bow, London E3 3PA

T: 020 7515 0333

E: dental@takara.co.uk

Manchester Showroom

7 Quays Reach, Carolina Way, Salford, Manchester M50 2ZY

T: 0161 743 9992

E: dental@takara.co.uk

All showroom visits by appointment only

Belmont

belmontdental.co.uk