



Ultrasoft
colours for Luxury &
Seamless upholstery

Advanced soft-to-touch fabric

 **Belmont**

Ultrasoft upholstery colours

for Belmont's Luxury and Seamless upholstery specifications




See product images in your preferred colours online at belmontdental.co.uk



Luxury Upholstery




Seamless Upholstery



291-7234 Emerald Isle




291-0001 White




291-5668 Silver Pearl




291-5666 Dove Grey




291-5040 Elephant




291-5708 Granite




291-5763 Charcoal



291-5806 Schooner



291-5607 Raven Wing



291-3666 Fudge




291-3205 Hot Chocolate




291-3779 Taupe




291-3602 Stone



291-1256 Papryus



291-3700 Ivory



291-3719 Champagne




291-3850 Sand



291-5250 Sassy



291-8243 Apricot




291-8233 Autumn Leaf



291-1347 Poppy




291-1176 Red



291-1312 Chianti




291-7214 Boysenberry



291-9385 Plum



291-2551 Baltic



291-2556 Admiral




291-2685 Bayou




291-2552 Riviera




291-2553 Cyan



291-2527 Lagoon



291-4351 Shetland



291-4314 Herbal Mist



291-4430 Eucalyptus



291-4454 Sonora



291-4460 Parrot



291-5200 Citrus

Ultrasoft

Soft-to-touch upholstery fabric
synonymous with luxury and durability

GENERAL SPECIFICATIONS

Composition:

100% Polycarbonate Polyurethane Surface
100% Reinforced Rayon Backing
containing TENCEL™ Modal fibers

Weight:

Approx. 9.7 oz. Per Square Yard
Approx. 330 Grams Per Square Metre

Thickness:

Approx. 0.95 Millimetres

Width:

54 Inches / 137 Centimetres

Repeat:

None

PERFORMANCE

Abrasion:

ASTM D4157 / 400,000 Double Rubs (Wyzenbeek)*
BS EN ISO 12947-2 / 180,000 Rubs (Martindale)*
*meets or exceeds

Hydrolysis:

14+ Weeks
ASTM D3690 Sect. 6.11 (Hydrolytic Stability)
ISO 1419 Method C (Accelerated Ageing)

Adhesion of Coating: ASTM D751
Tear Strength: ASTM D2261
Breaking Strength: ASTM D751
Antimicrobial: ASTM E2180
Tensile Strength: BS EN ISO 1421

COLOURFASTNESS

Light: AATCC 16.3, 200 Hours, Grade 4
ISO 105 B02, Grade 6

Crocking: AATCC 8, Grade 4-5 Wet & Dry

FLAMMABILITY REQUIREMENTS

- ASTM E84 (Adhered) – Class A
- CA TB 117-2013 – Pass
- FMVSS 302 & 95/28 Annex IV – Pass
- NFPA 260 – Class 1
- UFAC – Class 1

These fabrics, when treated and/or used with other compatible furniture components, have been shown to perform under the fire test conditions of:

- ASTM E84 (Unadhered)
- BS 5852 Source 0, 1 & 5 (Cig, Match & Crib 5)
- IMO FTPC Part 8
- BS EN 1021 Parts 1 & 2

For more information on certifications and standards, sustainability, as well as cleaning and care instructions, please visit: ultrafabricsinc.com

Made in Japan

 **Belmont**

Perfecting the art of dentistry



For any enquiry, do not hesitate to contact us directly on 020 7515 0333 or visit: belmontdental.co.uk

2025

Upholstery Disclaimer:

Belmont Luxury & Seamless Upholstery is supplied by a third party provider, Ultrasoft Fabrics. Whilst precautions are taken to ensure we can supply spare upholstered parts for your Belmont chair, Belmont cannot accept liability if a colour is discontinued at short notice by Ultrasoft Fabrics.